

1/8" fittings, 1.0 mm ports (.040")

MW Type

SPECS
200 psi gas
200°C max
 Nitronic 60 body
 Valcon E rotor

Includes 2" standoff.
 Ask about closemount
 assembly if valve will
 not be heated.

Standard electric actuators:
 110 VAC for USA
 110/230 VAC to 24 VDC power supply for international
 Microelectric actuators:
 24 VDC (includes a 110/230 VAC to 24 VDC power supply)

Low pressure

STF
Trap/ flow-throw

1/8" **1.0 mm**

OPTIONS

- 4 and 8 positions available
- 3", 4", and 6" standoffs
- Materials: Hastelloy C, Inconel 600, Monel 400, Nickel 200, Nitronic 50, Titanium (see pages 254-255)
- Larger bore available except 16 position

	6 Position		10 Position		12 Position		16 Position	
	Prod No	Price	Prod No	Price	Prod No	Price	Prod No	Price
Manual (not recommended)	2STF6MWE		2STF10MWE		2STF12MWE		2STF16MWE	
With air actuator	A2STF6MWE		A2STF10MWE		A2STF12MWE		A2STF16MWE	
With standard elec actuator	E2STF6MWE		E2STF10MWE		E2STF12MWE		E2STF16MWE	
With microelectric actuator	EMT2STF6MWE		EMT2STF10MWE		EMT2STF12MWE		EMT2STF16MWE	
Replacement valve	DSTF6MWE		DSTF10MWE		DSTF12MWE		DSTF16MWE	
Replacement rotor	SSASTF6MWE		SSASTF10MWE		SSASTF12MWE		SSASTF16MWE	

