

## NEW UHPLC Nanovolume Injectors with 360 µm fittings

### NEW 20,000 psi UHPLC Nanovolume valves 360 micron fittings, 100 micron bore (.004")

Model C72MU

**20,000 psi**  
**360µm 100 µm**

Includes stainless  
360 micron fittings.

Microelectric actuator:  
24 VDC, with 110/230 VAC to 24 VDC power supply.

#### SPECS

**20,000 psi liq**  
**50°C max**  
Stainless w/ inert coating  
stator  
Valcon E3 rotor



**4 Port**

Prod No

Price



**6 Port**

Prod No

Price

#### Coated stainless stator

Manual	C72MU-4694	C72MU-4696
With pneumatic actuator	C72MU-4694A	C72MU-4696A
With standard electric actuator	C72MU-4694E	C72MU-4696E
With microelectric actuator	C72MU-4694EH	C72MU-4696EH
Replacement valve	C72MU-4694D	C72MU-4696D
Replacement rotor	C72M-46R4	C72M-46R6
Replacement stator	C72M-4C94	C72M-4C96

#### OPTIONS

- 150 micron (.006") bore
- Internal sample injector (4 - 20 nl)
- 10,000 and 15,000 psi versions available



**Model C72MU**  
**360 micron fittings**  
(Model C72MX is similar)

### NEW 15,000 psi UHPLC Nanovolume valves 360 micron fittings, 150 micron bore (.006")

Model C72MX

**15,000 psi**  
**360µm 150 µm**

Includes stainless  
360 micron fittings.

Microelectric actuator:  
24 VDC, with 110/230 VAC to 24 VDC power supply.

#### SPECS

**15,000 psi liq**  
**50°C max**  
Stainless w/ inert coating  
stator  
Valcon E3 rotor



**6 Port**

Prod No

Price



**10 Port**

Prod No

Price

#### Coated stainless stator

Manual	C72MX-6696	C72MX-6690
With pneumatic actuator	C72MX-6696A	C72MX-6690A
With standard electric actuator	C72MX-6696E	C72MX-6690E
With microelectric actuator	C72MX-6696EH	C72MX-6690EH
Replacement valve	C72MX-6696D	C72MX-6690D
Replacement rotor	C72M-66R6	C72M-66R0
Replacement stator	C72M-6C96	C72M-6C90

#### OPTIONS

- 100 micron (.004") bore
- Internal sample injector (4 - 20 nl)
- 10,000 psi version available
- 4 and 8 port versions available

#### MORE INFORMATION

360 micron nanovolume fittings . . . . . page 58