








UHPLC Nanovolume Injectors with 1/32" Valco stainless fittings *NEW*

**NEW 15,000 psi UHPLC Nanovolume valves
1/32" Valco stainless fittings, 150 micron bore (.006")**

Model C72NX

<p>SPECS 15,000 psi liq 50°C max Stainless w/ inert coating stator Valcon E3 rotor</p>	<p>Includes 1/32" Valco stainless steel fittings.</p> 	<p>Microelectric actuator: 24 VDC, with 110/230 VAC to 24 VDC power supply.</p>	<p>* The 6 port valve includes a 5 µl loop of the stator material.</p>	15,000 psi
				1/32"

<p>OPTIONS</p> <ul style="list-style-type: none"> ■ 100 micron (.004") bore ■ 250 micron (.010") bore ■ 10,000 psi version available ■ 4 and 8 port versions available 	<p>Coated stainless stator</p> <p>Manual With pneumatic actuator</p> <p>With standard electric actuator With microelectric actuator</p> <p>Replacement valve Replacement rotor Replacement stator</p>	<table border="0"> <tr> <td></td> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> <tr> <td></td> <td style="text-align: center;">6 Port *</td> <td style="text-align: center;">10 Port</td> </tr> <tr> <td style="text-align: center;"><i>Prod No</i></td> <td style="text-align: center;"><i>Price</i></td> <td style="text-align: center;"><i>Prod No</i></td> <td style="text-align: center;"><i>Price</i></td> </tr> </table>					6 Port *	10 Port	<i>Prod No</i>	<i>Price</i>	<i>Prod No</i>	<i>Price</i>				
																
	6 Port *	10 Port														
<i>Prod No</i>	<i>Price</i>	<i>Prod No</i>	<i>Price</i>													
	<table border="0"> <tr> <td>C72NX-6696</td> <td>C72NX-6690</td> </tr> <tr> <td>C72NX-6696A</td> <td>C72NX-6690A</td> </tr> <tr> <td>C72NX-6696E</td> <td>C72NX-6690E</td> </tr> <tr> <td>C72NX-6696EH</td> <td>C72NX-6690ED</td> </tr> <tr> <td>C72NX-6696D</td> <td>C72NX-6690D</td> </tr> <tr> <td>C72N-66R6</td> <td>C72N-66R0</td> </tr> <tr> <td>C72N-6C96</td> <td>C72N-6C90</td> </tr> </table>	C72NX-6696	C72NX-6690	C72NX-6696A	C72NX-6690A	C72NX-6696E	C72NX-6690E	C72NX-6696EH	C72NX-6690ED	C72NX-6696D	C72NX-6690D	C72N-66R6	C72N-66R0	C72N-6C96	C72N-6C90	
C72NX-6696	C72NX-6690															
C72NX-6696A	C72NX-6690A															
C72NX-6696E	C72NX-6690E															
C72NX-6696EH	C72NX-6690ED															
C72NX-6696D	C72NX-6690D															
C72N-66R6	C72N-66R0															
C72N-6C96	C72N-6C90															



Model C72NX
1/32" Valco stainless fittings

Sample loops for C72NX valves

Each stainless loop includes two stainless 1/32" Cheminert nanovolume fittings.



Volume	Stainless	
	Prod No	Price
1 µl	CSLN1K	\$35.00
2 µl	CSLN2K	45.00
5 µl	CSLN5K	52.50



Model C74NX
1/32" Valco stainless fittings

**NEW 15,000 psi UHPLC Nanovolume internal sample injectors
1/32" Valco stainless fittings, 150 micron bore (.006")**

Model C74NX

<p>SPECS 15,000 psi liq 50°C max Stainless w/ inert coating stator Valcon E3 rotor</p>	<p>Includes 1/32" Valco stainless steel fittings.</p> 	<p>Microelectric actuator: 24 VDC, with 110/230 VAC to 24 VDC power supply.</p>		15,000 psi
				Internal sample

<p>OPTIONS</p> <ul style="list-style-type: none"> ■ 250 micron (.010") bore ■ 10,000 and 20,000 psi versions available 	<p>Coated stainless stator</p> <p>Manual With pneumatic actuator</p> <p>With standard electric actuator With microelectric actuator</p> <p>Replacement valve Replacement rotor Replacement stator</p>	<table border="0"> <tr> <td></td> <td style="text-align: center;">4 nanoliters</td> <td></td> <td style="text-align: center;">10 nanoliters</td> <td></td> <td style="text-align: center;">20 nanoliters</td> </tr> <tr> <td style="text-align: center;"><i>Prod No</i></td> <td style="text-align: center;"><i>Price</i></td> <td style="text-align: center;"><i>Prod No</i></td> <td style="text-align: center;"><i>Price</i></td> <td style="text-align: center;"><i>Prod No</i></td> <td style="text-align: center;"><i>Price</i></td> </tr> </table>		4 nanoliters		10 nanoliters		20 nanoliters	<i>Prod No</i>	<i>Price</i>	<i>Prod No</i>	<i>Price</i>	<i>Prod No</i>	<i>Price</i>									
			4 nanoliters		10 nanoliters		20 nanoliters																
<i>Prod No</i>	<i>Price</i>	<i>Prod No</i>	<i>Price</i>	<i>Prod No</i>	<i>Price</i>																		
	<table border="0"> <tr> <td>C74NX-6694-.004</td> <td>C74NX-6694-.01</td> <td>C74NX-6694-.02</td> </tr> <tr> <td>C74NX-6694-.004A</td> <td>C74NX-6694-.01A</td> <td>C74NX-6694-.02A</td> </tr> <tr> <td>C74NX-6694-.004E</td> <td>C74NX-6694-.01E</td> <td>C74NX-6694-.02E</td> </tr> <tr> <td>C74NX-6694-.004EH</td> <td>C74NX-6694-.01EH</td> <td>C74NX-6694-.02EH</td> </tr> <tr> <td>C74NX-6694-.004D</td> <td>C74NX-6694-.01D</td> <td>C74NX-6694-.02D</td> </tr> <tr> <td>C74N-66R-.004</td> <td>C74N-66R-.01</td> <td>C74N-66R-.02</td> </tr> <tr> <td>C74N-6C9</td> <td>C74N-6C9</td> <td>C74N-6C9</td> </tr> </table>	C74NX-6694-.004	C74NX-6694-.01	C74NX-6694-.02	C74NX-6694-.004A	C74NX-6694-.01A	C74NX-6694-.02A	C74NX-6694-.004E	C74NX-6694-.01E	C74NX-6694-.02E	C74NX-6694-.004EH	C74NX-6694-.01EH	C74NX-6694-.02EH	C74NX-6694-.004D	C74NX-6694-.01D	C74NX-6694-.02D	C74N-66R-.004	C74N-66R-.01	C74N-66R-.02	C74N-6C9	C74N-6C9	C74N-6C9	
C74NX-6694-.004	C74NX-6694-.01	C74NX-6694-.02																					
C74NX-6694-.004A	C74NX-6694-.01A	C74NX-6694-.02A																					
C74NX-6694-.004E	C74NX-6694-.01E	C74NX-6694-.02E																					
C74NX-6694-.004EH	C74NX-6694-.01EH	C74NX-6694-.02EH																					
C74NX-6694-.004D	C74NX-6694-.01D	C74NX-6694-.02D																					
C74N-66R-.004	C74N-66R-.01	C74N-66R-.02																					
C74N-6C9	C74N-6C9	C74N-6C9																					

MORE INFORMATION
1/32" Valco fittings pages 10, 12
1/16" nanovolume injectors
C72X156
C74X157

CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034