

Nanovolume Injectors with 1/32" Cheminert Fittings

5,000 psi Nanovolume internal sample injector, 1/32" Cheminert fittings, 100 micron ports (.004")

Model CN4

SPECS
5000 psi liq
50°C max
PAEK stator
Valcon E rotor

Includes 1/32" PEEK Cheminert nanovolume fittings.



Microelectric actuator:
24 VDC, with 110/230 VAC to 24 VDC power supply.



5,000 psi
Internal sample
1/32" **100 µm**

OPTIONS

- 150 micron (.006") bore

Sample volume	4 nanoliters		10 nanoliters		20 nanoliters	
	Prod No	Price	Prod No	Price	Prod No	Price
PAEK stator						
Manual	CN4-4344-.004		CN4-4344-.01		CN4-4344-.02	
With microelectric actuator	CN4-4344-.004EH		CN4-4344-.01EH		CN4-4344-.02EH	
Replacement valve	CN4-4344-.004D		CN4-4344-.01D		CN4-4344-.02D	
Replacement rotor	CN4-43R-.004		CN4-43R-.01		CN4-43R-.02	
Replacement stator	CN4-4C4I		CN4-4C4I		CN4-4C4I	



Model CN4
1/32" Cheminert fittings

MORE INFORMATION

1/32" PEEK Cheminert fitting (nut with collapsible ferrule) p.59

CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034