

Microbore HPLC

Microbore valves, 1/16" Valco fittings, 0.25 mm ports (.010")

Model C2

5,000 psi

Microbore

1/16" 0.25 mm

Includes stainless steel nuts and ferrules of the stator material. Valves with PEEK stators have PEEK nuts and ferrules.

Standard electric actuator:
110 VAC for USA
110/230 VAC to 24 VDC power supply for international.

Microelectric actuator:
24 VDC, with 110/230 VAC to 24 VDC power supply.

* The 6 port valve includes a 5 µl loop of the stator material.

SPECS
5000 psi liq
75°C max
Metal stator
Valcon H rotor

5000 psi liq
50°C max
PAEK stator
Valcon E rotor



4 Port

Prod No Price



6 Port*

Prod No Price



8 Port

Prod No Price



10 Port

Prod No Price

N60 stainless stator

Manual	C2-1004	C2-1006	C2H-1008	C2H-1000
With pneumatic actuator	C2-1004A	C2-1006A	C2H-1008A	C2H-1000A
With standard electric actuator	C2-1004E	C2-1006E	C2H-1008E	C2H-1000E
With microelectric actuator	C2-1004EH	C2-1006EH	C2H-1008EH	C2H-1000EH
Replacement valve	C2-1004D	C2-1006D	C2H-1008D	C2H-1000D
Replacement rotor	C2-10R4	C2-10R6	C2-10R8H	C2-10R0H
Replacement stator	C-1C04	C-1C06	C-1C08H	C-1C00H

PAEK stator

Manual	C2-1344	C2-1346	C2H-1348	C2H-1340
With pneumatic actuator	C2-1344A	C2-1346A	C2H-1348A	C2H-1340A
With standard electric actuator	C2-1344E	C2-1346E	C2H-1348E	C2H-1340E
With microelectric actuator	C2-1344EH	C2-1346EH	C2H-1348EH	C2H-1340EH
Replacement valve	C2-1344D	C2-1346D	C2H-1348D	C2H-1340D
Replacement rotor	C2-13R4	C2-13R6	C2-13R8H	C2-13R0H
Replacement stator	C-1C44	C-1C46	C-1C48H	C-1C40H

Titanium stator

Manual	C2-1034	C2-1036	C2H-1038	C2H-1030
With pneumatic actuator	C2-1034A	C2-1036A	C2H-1038A	C2H-1030A
With standard electric actuator	C2-1034E	C2-1036E	C2H-1038E	C2H-1030E
With microelectric actuator	C2-1034EH	C2-1036EH	C2H-1038EH	C2H-1030EH
Replacement valve	C2-1034D	C2-1036D	C2H-1038D	C2H-1030D
Replacement rotor	C2-10R4	C2-10R6	C2-10R8H	C2-10R0H
Replacement stator	C-1C34	C-1C36	C-1C38H	C-1C30H

OPTIONS

- Continuous flow version is available as Model C6. See page 160.
- Hastelloy C stators
- Loop fill port assembly for injection from front of the valve. See page 41.
- 0.15 mm (0.006") bore



Order loops from page 159.



Model C2
1/16" ZDV fittings

MORE INFORMATION

- Actuators
- Air page 195
 - Microelectric .. 188-189
 - Standard electric ... 193
- Materials
- Metals..... 254-255
 - Polymers 256
 - Valve rotors..... 257
- Nuts
- Metal..... 10
 - PEEK 63
- Ferrules
- Metal..... 12
 - PEEK 63
- Standoff assemblies 205-207