

**Analytical valves,
1/16" Valco fittings, 0.40 mm ports (.016")**

Model C2

SPES
5000 psi liq
75°C max
 Metal stator
 Valcon H rotor

5000 psi liq
50°C max
 PAEK stator
 Valcon E rotor

Includes stainless steel nuts and ferrules of the stator material. Valves with PAEK stators have PEEK nuts and ferrules.

Standard electric actuator:
 110 VAC for USA
 110/230 VAC to 24 VDC power supply for international.
 Microelectric actuator:
 24 VDC, with 110/230 VAC to 24 VDC power supply.

* The 6 port valve includes a 20 µl loop of the stator material.

5,000 psi

Analytical

1/16"

0.40 mm

OPTIONS

- Continuous flow version is available as Model C6. See page 163.
- Hastelloy C stators
- Semi-prep version with 0.75 mm ports (.030") available
- Loop fill port assembly for injection from front of the valve. See page 41.



Order loops from page 159.

OPTIONAL FLOWPATH

Model C2 6 port valves can also be ordered with a dual 3-way rotor, as described in EPA Method 555.

To specify this flowpath, substitute "6X" for "6" in the valve or rotor product number.



AUTOSAMPLER REPLACEMENT VALVES

The Cheminert Model C2 6 port valve is an excellent replacement for the valve originally supplied in many autosamplers, including autosamplers manufactured by Beckman, Gilson, Spark-Holland, CTC, Thermo Fisher, and Varian.

Call technical support to determine which replacement is best for your application.



4 Port



6 Port*



8 Port



10 Port

Prod No Price Prod No Price Prod No Price Prod No Price

N60 stainless stator

Manual	C2-2004	C2-2006	C2H-2008	C2H-2000
With pneumatic actuator	C2-2004A	C2-2006A	C2H-2008A	C2H-2000A
With standard electric actuator	C2-2004E	C2-2006E	C2H-2008E	C2H-2000E
With microelectric actuator	C2-2004EH	C2-2006EH	C2H-2008EH	C2H-2000EH
Replacement valve	C2-2004D	C2-2006D	C2H-2008D	C2H-2000D
Replacement rotor	C2-20R4	C2-20R6	C2-20R8H	C2-20R0H
Replacement stator	C-2C04	C-2C06	C-2C08H	C-2C00H

Replacement valve	C2-2344D	C2-2346D	C2H-2348D	C2H-2340D
Replacement rotor	C2-23R4	C2-23R6	C2-23R8H	C2-23R0H
Replacement stator	C-2C44	C-2C46	C-2C48H	C-2C40H

Titanium stator

Manual	C2-2034	C2-2036	C2H-2038	C2H-2030
With pneumatic actuator	C2-2034A	C2-2036A	C2H-2038A	C2H-2030A
With standard electric actuator	C2-2034E	C2-2036E	C2H-2038E	C2H-2030E
With microelectric actuator	C2-2034EH	C2-2036EH	C2H-2038EH	C2H-2030EH
Replacement valve	C2-2034D	C2-2036D	C2H-2038D	C2H-2030D
Replacement rotor	C2-20R4	C2-20R6	C2-20R8H	C2-20R0H
Replacement stator	C-2C34	C-2C36	C-2C38H	C-2C30H



**Model C2
1/16" ZDV fittings**