

Analytical HPLC

Analytical internal sample injector, 1/16" Valco fittings, 0.25 mm ports (.010")

Model C4

5,000 psi

Analytical

Internal sample

1/16" 0.25 mm



Includes stainless steel nuts and ferrules of the stator material. Valves with PAEK stators have PEEK nuts and ferrules.

Standard electric actuator:
110 VAC for USA
110/230 VAC to 24 VDC power supply for international.

Microelectric actuator:
24 VDC, with 110/230 VAC to 24 VDC power supply.

SPECS
5000 psi liq
75°C max
Metal stator
Valcon H rotor

5000 psi liq
50°C max
PAEK stator
Valcon E rotor

Sample volume	0.1 µl		0.2 µl		0.5 µl	
	Prod No	Price	Prod No	Price	Prod No	Price
N60 stainless stator						
Manual	C4-1004-.1		C4-1004-.2		C4-1004-.5	
With pneumatic actuator	C4-1004-.1A		C4-1004-.2A		C4-1004-.5A	
With standard electric actuator	C4-1004-.1E		C4-1004-.2E		C4-1004-.5E	
With microelectric actuator	C4-1004-.1EH		C4-1004-.2EH		C4-1004-.5EH	
Replacement valve	C4-1004-.1D		C4-1004-.2D		C4-1004-.5D	
Replacement rotor	C4-10R-.1		C4-10R-.2		C4-10R-.5	
Replacement stator	C4-1C0		C4-1C0		C4-1C0	
PAEK stator						
Manual	C4-1344-.1		C4-1344-.2		C4-1344-.5	
With pneumatic actuator	C4-1344-.1A		C4-1344-.2A		C4-1344-.5A	
With standard electric actuator	C4-1344-.1E		C4-1344-.2E		C4-1344-.5E	
With microelectric actuator	C4-1344-.1EH		C4-1344-.2EH		C4-1344-.5EH	
Replacement valve	C4-1344-.1D		C4-1344-.2D		C4-1344-.5D	
Replacement rotor	C4-13R-.1		C4-13R-.2		C4-13R-.5	
Replacement stator	C4-1C4		C4-1C4		C4-1C4	
Titanium stator						
Manual	C4-1034-.1		C4-1034-.2		C4-1034-.5	
With pneumatic actuator	C4-1034-.1A		C4-1034-.2A		C4-1034-.5A	
With standard electric actuator	C4-1034-.1E		C4-1034-.2E		C4-1034-.5E	
With microelectric actuator	C4-1034-.1EH		C4-1034-.2EH		C4-1034-.5EH	
Replacement valve	C4-1034-.1D		C4-1034-.2D		C4-1034-.5D	
Replacement rotor	C4-10R-.1		C4-10R-.2		C4-10R-.5	
Replacement stator	C4-1C3		C4-1C3		C4-1C3	

OPTIONS

- .05 µl sample volumes are also available.
- Loop fill port assembly for injection from front of the valve.
See page 41.



Model C4
1/16" ZDV fittings

MORE INFORMATION

- Actuators
- Air page 195
 - Microelectric .. 188-189
 - Standard electric ... 193
- Materials
- Metals..... 254-255
 - Polymers 256
 - Valve rotors..... 257
- Standoff assemblies 205-207