





Valves with 1/4-28 fitting details for 1/16" tubing, 0.75 mm ports (.030")

Model C22

SPECS 100 psi gas/ 250 psi liq 75°C max PPS stator Valcon E2 rotor	Includes multicolored Cheminert 1/4-28 flangeless fittings for 1/16" tubing.	Standard electric actuator: 110 VAC for USA 110/230 VAC to 24 VDC power supply for international. Microelectric actuator: 24 VDC, with 110/230 VAC to 24 VDC power supply.	Sample loops are not included with valves. Order separately.	Low pressure			
				1/4-28 Internal			
				1/16"	0.75 mm		
							
				4 Port	6 Port	8 Port	10 Port
				<i>Prod No</i> <i>Price</i>	<i>Prod No</i> <i>Price</i>	<i>Prod No</i> <i>Price</i>	<i>Prod No</i> <i>Price</i>
Manual				C22-3184	C22-3186	C22-3188	C22-3180
With pneumatic actuator				C22-3184A	C22-3186A	C22-3188A	C22-3180A
With standard electric actuator				C22-3184E	C22-3186E	C22-3188E	C22-3180E
With microelectric actuator				C22-3184EH	C22-3186EH	C22-3188EH	C22-3180EH
Replacement valve				C22-3184D	C22-3186D	C22-3188D	C22-3180D
Replacement rotor				C22-314	C22-316	C22-318	C22-310
Replacement stator				C22-384	C22-386	C22-388	C22-380

Valves with 1/4-28 fitting details for 1/8" tubing, 1.50 mm ports (.060")

Model

SPECS 100 psi gas/ 250 psi liq 75°C max PPS stator Valcon E2 rotor	Includes multicolored Cheminert 1/4-28 flangeless fittings for 1/8" tubing.	Standard electric actuator: 110 VAC for USA 110/230 VAC to 24 VDC power supply for international. Microelectric actuator: 24 VDC, with 110/230 VAC to 24 VDC power supply.	Sample loops are not included with valves. Order separately.	Low pressure			
				1/4-28 Internal			
				1/8"	1.50 mm		
				4 Port	6 Port	8 Port	10 Port
				<i>Prod No</i> <i>Price</i>	<i>Prod No</i> <i>Price</i>	<i>Prod No</i> <i>Price</i>	<i>Prod No</i> <i>Price</i>
Manual				C22-6184	C22-6186	C22-6188	C22-6180
With pneumatic actuator				C22-6184A	C22-6186A	C22-6188A	C22-6180A
With standard electric actuator				C22-6184E	C22-6186E	C22-6188E	C22-6180E
With microelectric actuator				C22-6184EH	C22-6186EH	C22-6188EH	C22-6180EH
Replacement valve				C22-6184D	C22-6186D	C22-6188D	C22-6180D
Replacement rotor				C22-614	C22-616	C22-618	C22-610
Replacement stator				C22-684	C22-686	C22-688	C22-680

Sample loops

for Model C22

Loops include flangeless fittings with white color nuts. Loops smaller than 500 µl are made from 1/16" OD tubing; loops 500 µl or bigger are made from 1/8" OD tubing.

Volume	FEP		PTFE		PEEK	
	Prod No	Price	Prod No	Price	Prod No	Price
20 µl	CFSL20FEP		CFSL20TF	\$17.50	CFSL20PK	
50 µl	CFSL50FEP		CFSL50TF		CFSL50PK	
100 µl	CFSL100FEP		CFSL100TF		CFSL100PK	
250 µl	CFSL250FEP		CFSL250TF		CFSL250PK	
500 µl	CFSL500FEP		CFSL500TF		CFSL500PK	
1 ml	CFSL1KFEP		CFSL1KTF		CFSL1KPK	
2 ml	CFSL2KFEP		CFSL2KTF		CFSL2KPK	



Model C22
1/4-28 fittings