

**Internal sample injectors, 1/4-28 for 1/16" tubing, 0.50 mm ports (.020")**

*Model C24*

**SPECS**  
**100 psi gas/ 250 psi liq**  
**75°C max**  
 PPS stator  
 Valcon E2 rotor

Includes multicolored Cheminert 1/4-28 flangeless fittings for 1/16" tubing.

Standard electric actuator:  
 110 VAC for USA  
 110/230 VAC to 24 VDC power supply for international.  
 Microelectric actuator:  
 24 VDC, with 110/230 VAC to 24 VDC power supply.



- Low pressure**
- Internal sample**
- 1/4-28 Internal**
- 1/16" 0.50 mm**

**OPTIONS**

- 0.2 µl sample volumes are also available.
- Purge option
- Other polymeric rotors and stators are available. Consult the factory for prices and information.

Sample volume	0.5 µl		1 µl		2 µl	
	Prod No	Price	Prod No	Price	Prod No	Price
Manual	C24-2184-.5		C24-2184-1		C24-2184-2	
With pneumatic actuator	C24-2184-.5A		C24-2184-1A		C24-2184-2A	
With standard electric actuator	C24-2184-.5E		C24-2184-1E		C24-2184-2E	
With microelectric actuator	C24-2184-.5EH		C24-2184-1EH		C24-2184-2EH	
Replacement valve	C24-2184-.5D		C24-2184-1D		C24-2184-2D	
Replacement rotor	C24-10R-.5		C24-10R-1		C24-10R-2	
Replacement stator	C24-1C8		C24-1C8		C24-1C8	



**Model C24**  
 1/4-28 fittings