

NEW Selectors – Nanovolume UHPLC

NEW 15,000 psi UHPLC Nanovolume selectors, 1/32" Valco fittings, 150 micron ports (.006")


Model C75NX

- 15,000 psi**
- Stream selector**
- 1/32"**
- 150 µm**

Includes 1/32" Valco stainless steel fittings.



Manual version not available.
 Standard electric actuator:
 110 VAC for USA
 110/230 VAC to 24 VDC power supply for international.
 Microelectric actuator:
 24 VDC, with 110/230 VAC to 24 VDC power supply.



SPECS

15,000 psi liq
50°C max
 Stainless stator with inert coating
 Valcon E3 rotor

OPTIONS

- 100 micron (.004") bore
- 250 micron (.010") bore
- 10,000 and 20,000 psi versions available
- 4 positions

	6 Position		8 Position		10 Position	
	Prod No	Price	Prod No	Price	Prod No	Price
Coated stainless stator						
With pneumatic actuator	C75NX-6696A		C75NX-6698A		C75NX-6690A	
With standard electric actuator	C75NX-6696E		C75NX-6698E		C75NX-6690E	
With microelectric actuator	C75NX-6696EMH		C75NX-6698EMT		C75NX-6690EMT	
Replacement valve	C75NX-6696D		C75NX-6698D		C75NX-6690D	
Replacement rotor	C75N-66R6		C75N-66R8		C75N-66R0	
Replacement stator	C75N-6C96		C75N-6C98		C75N-6C90	



Model C75NX
 1/32" Valco stainless fittings