

**NEW 10,000 psi UHPLC microbore selectors,
1/16" Valco fittings, 250 micron ports (.010")**

Model C75H

SPECS
10,000 psi liq
50°C max
 Stainless stator
 with inert coating
 Valcon E3 rotor

Includes 1/16" Valco
 stainless steel fittings.



Manual version not available.
 Standard electric actuator:
 110 VAC for USA
 110/230 VAC to 24 VDC power supply for
 international.
 Microelectric actuator:
 24 VDC, with 110/230 VAC to 24 VDC
 power supply.

10,000 psi

Stream selector

1/16"

250 µm

OPTIONS

- 150 micron (.006") bore
- 15,000 psi version available
- 4 positions

	6 Position		8 Position		10 Position	
	Prod No	Price	Prod No	Price	Prod No	Price
Coated stainless stator						
With pneumatic actuator	C75H-1696A		C75H-1698A		C75H-1690A	
With standard electric actuator	C75H-1696E		C75H-1698E		C75H-1690E	
With microelectric actuator	C75H-1696EMH		C75H-1698EMT		C75H-1690EMT	
Replacement valve	C75H-1696D		C75H-1698D		C75H-1690D	
Replacement rotor	C75-16R6		C75-16R8		C75-16R0	
Replacement stator	C75-1C96		C75-1C98		C75-1C90	



Model C75H
 1/16" Valco stainless
 fittings