


### Stream selector, 1/4-28 fittings for 1/16" tubing, 0.75 mm ports (.030")

Model C25

**SPECS**  
**100 psi gas/ 250 psi liq**  
**75°C max**  
 PPS stator  
 Valcon E2 rotor

Includes multicolored Cheminert 1/4-28 flangeless fittings for 1/16" tubing.

Standard electric actuator:  
 110 VAC for USA  
 110/230 VAC to 24 VDC power supply for international.  
 Microelectric actuator:  
 24 VDC, with 110/230 VAC to 24 VDC power supply.



**Low pressure**  
**Stream selector**  
**1/4-28 Internal**  
**1/16" 0.75 mm**

**OPTIONS**

- 2", 3", 4", and 6" standoffs
- CTFE stator

	4 Position		6 Position		8 Position		10 Position	
	Prod No	Price	Prod No	Price	Prod No	Price	Prod No	Price
Manual	C25-3184		C25-3186		C25-3188		C25-3180	
With pneumatic act.	C25-3184A		C25-3186A		C25-3188A		C25-3180A	
With std electric act.	C25-3184E		C25-3186E		C25-3188E		C25-3180E	
With microelec act.	C25-3184EMH		C25-3186EMH		C25-3188EMH		C25-3180EMH	
Replacement valve	C25-3184D		C25-3186D		C25-3188D		C25-3180D	
Replacement rotor	C25-314		C25-316		C25-318		C25-310	
Replacement stator	C25-384		C25-386		C25-388		C25-380	


### Stream selector, 1/4-28 fittings for 1/8" tubing, 1.50 mm ports (.060")

Model C25

**SPECS**  
**100 psi gas/ 250 psi liq**  
**75°C max**  
 PPS stator  
 Valcon E2 rotor

Includes multicolored Cheminert 1/4-28 flangeless fittings for 1/8" tubing.

Standard electric actuator:  
 110 VAC for USA  
 110/230 VAC to 24 VDC power supply for international.  
 Microelectric actuator:  
 24 VDC, with 110/230 VAC to 24 VDC power supply.



**Low pressure**  
**Stream selector**  
**1/4-28 Internal**  
**1/8" 1.50 mm**

**OPTIONS**

- 2", 3", 4", and 6" standoffs
- CTFE stator

**MORE INFORMATION**

Actuators  
 Air ..... page 194  
 Microelectric .. 190-191  
 Standard electric ... 193  
 Materials  
 Metals..... 254-255  
 Polymers ..... 256  
 Valve rotors..... 257  
 Standoff assemblies .... 205-207

**OPTIONAL FLOWPATH**

**Model C25F** is the flow-through version of C25. (See discussion on facing page.) 3, 4, 5, 6, and 7 positions are available.

Consult the factory for C25F prices and information.

	4 Position		6 Position		8 Position		10 Position	
	Prod No	Price	Prod No	Price	Prod No	Price	Prod No	Price
Manual	C25-6184		C25-6186		C25-6188		C25-6180	
With pneumatic act.	C25-6184A		C25-6186A		C25-6188A		C25-6180A	
With std electric act.	C25-6184E		C25-6186E		C25-6188E		C25-6180E	
With microelec act.	C25-6184EMH		C25-6186EMH		C25-6188EMH		C25-6180EMH	
Replacement valve	C25-6184D		C25-6186D		C25-6188D		C25-6180D	
Replacement rotor	C25-614		C25-616		C25-618		C25-610	
Replacement stator	C25-684		C25-686		C25-688		C25-680	



**Model C25**  
**10 position**  
**1/4-28 fittings**