


Stream selector, 1/2-20 fittings for 1/4" tubing, 4.6 mm ports (.180")

Model C45

SPECS
100 psi liq
50°C max
 PPS stator
 Valcon TF rotor

Manual version not available.
 Includes Cheminert 1/2-20 flangeless fittings for 1/4" tubing, Delrin nuts and CTFE ferrules.

Standard electric actuator:
 110 VAC for USA
 110/230 VAC to 24 VDC power supply for international.
 Microelectric actuator:
 24 VDC, with 110/230 VAC to 24 VDC power supply.



Low pressure

Stream selector

1/2-20 Internal

1/4" 4.6 mm

OPTIONS

- 2", 3", 4", and 6" standoffs
- Consult the factory for optional materials.
- 8 position selectors are available with 3 mm (.120") ports


	4 Position		6 Position	
	Prod No	Price	Prod No	Price
With pneumatic actuator	C45-9784A		C45-9786A	
With std electric actuator	C45-9784E		C45-9786E	
With microelectric actuator	C45-9784EMT		C45-9786EMT	
Replacement valve	C45-9784D		C45-9786D	
Replacement rotor	C45-97R4		C45-97R6	

Fittings for C45 valves



For additional 1/2-20 fittings and adapters, see page 72.

	Prod No	Price
Delrin nut	CFL-4D	
CTFE nut	CFL-4KF	
PPS nut	CFL-4PPS	
CTFE ferrule	CFL-CB4KF-S	



Model C45
6 positions
1/2-20 fittings

MORE INFORMATION

- Actuators
 Air page 194
 Microelectric .. 190-191
 Standard electric ... 193
- Materials
 Metals..... 254-255
 Polymers 256
 Valve rotors..... 257
- Standoff assemblies 205-207