

Analytical centered port injector, 1/16" Valco fittings, 0.40 mm ports (.016")

Model C3

SPECS
5000 psi liq
75°C max
 Metal stator
 Valcon H rotor

5000 psi liq
50°C max
 PAEK stator
 Valcon E rotor

Includes stainless steel nuts and ferrules. Valves with PAEK stators have PEEK nuts and ferrules.

Includes one 20 µl loop of the stator material. Includes syringe fill port for 22 gauge 3/4" and 2" needle.


Standard electric actuator:
 110 VAC for USA
 110/230 VAC to 24 VDC power supply for international.

Microelectric actuator:
 24 VDC, with 110/230 VAC to 24 VDC power supply.




- 5,000 psi**
- Analytical**
- Centered port**
- 1/16"**
- 0.40 mm**

N60 stainless stator		PAEK stator	
Prod No	Price	Prod No	Price
Manual	C3-2006	C3-2346	
With pneumatic actuator	C3-2006A	C3-2346A	
With standard electric actuator	C3-2006E	C3-2346E	
With microelectric actuator	C3-2006EH	C3-2346EH	
Replacement valve	C3-2006D	C3-2346D	
Replacement rotor	C2-20R6	C2-23R6	
Replacement stator	C3-2C06	C3-2C46	



OPTIONS
 ■ Titanium and Hastelloy stators available



Analytical vertical port injector, 1/16" Valco fittings, 0.40 mm ports (.016")

Model C2V

SPECS
5000 psi liq
75°C max
 Metal stator
 Valcon H rotor

5000 psi liq
50°C max
 PAEK stator
 Valcon E rotor

Includes stainless steel nuts and ferrules. Valves with PAEK stators have PEEK nuts and ferrules.

Includes one 20 µl loop of the stator material.


Standard electric actuator:
 110 VAC for USA
 110/230 VAC to 24 VDC power supply for international.

Microelectric actuator:
 24 VDC, with 110/230 VAC to 24 VDC power supply.



- 5,000 psi**
- Analytical**
- Vertical port**
- 1/16"**
- 0.40 mm**

N60 stainless stator		PAEK stator	
Prod No	Price	Prod No	Price
Manual	C2V-2006	C2V-2346	
With pneumatic actuator	C2V-2006A	C2V-2346A	
With standard electric actuator	C2V-2006E	C2V-2346E	
With microelectric actuator	C2V-2006EH	C2V-2346EH	
Replacement valve	C2V-2006D	C2V-2346D	
Replacement rotor	C2-20R6	C2-23R6	
Replacement stator	C2V-2C06	C2V-2C46	



OPTIONS
 ■ Titanium and Hastelloy stators available

