

OEM – Low Pressure

NEW Integrated motor/valve, 1/16" Valco ZDV fittings, 0.75 mm ports (.030")

Model C62Z

- Low pressure
- Integrated
- 10-32 ZDV
- 1/16"
- 0.75 mm

Includes Valco ZDV PEEK nuts and ferrules.

Sample loops are not included with valves. Order separately.

SPECS

100 psi gas/ 250 psi liq
50°C max
PPS stator
Valcon E2 rotor

CE



4 Port

Prod No Price



6 Port

Prod No Price



8 Port

Prod No Price



10 Port

Prod No Price

With motor, sensor, & controller
With motor and sensor only

C62Z-3184I
C62Z-3184I-S

C62Z-3186I
C62Z-3186I-S


C62Z-3188I
C62Z-3188I-S

C62Z-3180I
C62Z-3180I-S

OPTIONS

- Other polymeric rotors and stators are available

Consult the factory for prices and information.



Model C62Z
1/16" ZDV fittings

Sample loops

for Model C62Z

Loops include PEEK nuts and ferrules. Loops less than 500 µl are made from 1/16" OD tubing; loops 500 µl or greater are made from 1/8" OD tubing with polymeric unions and 1/16" ends.



Volume	FEP		PTFE		PEEK	
	Prod No	Price	Prod No	Price	Prod No	Price
5 µl	CZSL5FEP		CZSL5TF		CZSL5PK	
10 µl	CZSL10FEP		CZSL10TF		CZSL10PK	
20 µl	CZSL20FEP		CZSL20TF		CZSL20PK	
50 µl	CZSL50FEP		CZSL50TF		CZSL50PK	
100 µl	CZSL100FEP		CZSL100TF		CZSL100PK	
250 µl	CZSL250FEP		CZSL250TF		CZSL250PK	
500 µl	CZSL500FEP		CZSL500TF		CZSL500PK	
1 ml	CZSL1KFEP		CZSL1KTF		CZSL1KPK	
2 ml	CZSL2KFEP		CZSL2KTF		CZSL2KPK	

MORE INFORMATION

Materials
Metals page 254-255
Polymers 256
Valve rotors 257

CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034