

OEM – Selectors – High Pressure


NEW Integrated motor/stream selector, 1/16" Valco ZDV fittings, 0.40 mm ports (.016")

Model C55

- 5,000 psi
- Integrated
- Stream selector
- 10-32 ZDV
- 1/16" 0.40 mm

Includes stainless steel nuts and ferrules of the stator material.

Valves with PAEK stators have PEEK nuts and ferrules.





SPECS
5000 psi liq
50°C max
 Metal stator
 Valcon H rotor

5000 psi liq
50°C max
 PAEK stator
 Valcon E rotor



	4 Position <i>Prod No</i>	6 Position <i>Prod No</i>	8 Position <i>Prod No</i>	10 Position <i>Prod No</i>
N60 stainless stator				
With integrated actuator	C55-2004I	C55-2006I	C55-2008I	C55-2000I
With motor/sensor only	C55-2004I-S	C55-2006I-S	C55-2008I-S	C55-2000I-S
With motor only	C55-2004IX	C55-2006IX	C55-2008IX	C55-2000IX
Replacement rotor	C5-20R4	C5-20R6	C5-20R8H	C5-20R0H
Replacement stator	C55-2C04	C55-2C06	C55-2C08	C55-2C00
PAEK stator				
With integrated actuator	C55-2344I	C55-2346I	C55-2348I	C55-2340I
With motor/sensor only	C55-2344I-S	C55-2346I-S	C55-2348I-S	C55-2340I-S
With motor only	C55-2344IX	C55-2346IX	C55-2348IX	C55-2340IX
Replacement rotor	C5-23R4	C5-23R6	C5-23R8H	C5-23R0H
Replacement stator	C55-2C44	C55-2C46	C55-2C48	C55-2C40

- OPTIONS**
- Optional bore:
0.25 mm (.010")
0.75 mm (.030")
 - 4 and 8 positions available



CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034