


### NEW Integrated motor/stream selector, 1/16" Valco ZDV fittings, 0.75 mm ports (.030")

Model C65Z

**SPECS**  
100 psi gas/ 250 psi liq  
50°C max  
PPS stator  
Valcon E2 rotor

Includes Valco ZDV  
PEEK nuts and ferrules.



**Low pressure**  
**Integrated**  
**Stream selector**  
**10-32 ZDV**  
**1/16" 0.75 mm**  
CE


|                            | <b>4 Position</b><br><i>Prod No</i> | <b>6 Position</b><br><i>Prod No</i> | <b>8 Position</b><br><i>Prod No</i> | <b>10 Position</b><br><i>Prod No</i> |
|----------------------------|-------------------------------------|-------------------------------------|-------------------------------------|--------------------------------------|
| With integrated actuator   | C65Z-3184I                          | C65Z-3186I                          | C65Z-3188I                          | C65Z-3180I                           |
| With motor and sensor only | C65Z-3184I-S                        | C65Z-3186I-S                        | C65Z-3188I-S                        | C65Z-3180I-S                         |

### NEW Integrated motor/stream stream selector, 1/4-28 fittings for 1/16" tubing, 0.75 mm ports (.030")

Model C65

**SPECS**  
100 psi gas/ 250 psi liq  
50°C max  
PPS stator  
Valcon E2 rotor

Includes multicolored  
Cheminert flangeless  
fittings for 1/16"  
tubing. See photo on facing  
page.



**Low pressure**  
**Integrated**  
**Stream selector**  
**1/4-28 Internal**  
**1/16" 0.75 mm**  
CE


|                            | <b>4 Position</b><br><i>Prod No</i> | <b>6 Position</b><br><i>Prod No</i> | <b>8 Position</b><br><i>Prod No</i> | <b>10 Position</b><br><i>Prod No</i> |
|----------------------------|-------------------------------------|-------------------------------------|-------------------------------------|--------------------------------------|
| With integrated actuator   | C65-3184I                           | C65-3186I                           | C65-3188I                           | C65-3180I                            |
| With motor and sensor only | C65-3184I-S                         | C65-3186I-S                         | C65-3188I-S                         | C65-3180I-S                          |

### NEW Integrated motor/stream stream selector, 1/4-28 fittings for 1/8" tubing, 1.50 mm ports (.060")

Model C65

**SPECS**  
100 psi gas/ 250 psi liq  
50°C max  
PPS stator  
Valcon E2 rotor

Includes multicolored  
Cheminert flangeless  
fittings for 1/8"  
tubing. See photo on facing  
page.



**Low pressure**  
**Integrated**  
**Stream selector**  
**1/4-28 Internal**  
**1/8" 1.50 mm**  
CE

|                            | <b>4 Position</b><br><i>Prod No</i> | <b>6 Position</b><br><i>Prod No</i> | <b>8 Position</b><br><i>Prod No</i> | <b>10 Position</b><br><i>Prod No</i> |
|----------------------------|-------------------------------------|-------------------------------------|-------------------------------------|--------------------------------------|
| With integrated actuator   | C65-6184I                           | C65-6186I                           | C65-6188I                           | C65-6180I                            |
| With motor and sensor only | C65-6184I-S                         | C65-6186I-S                         | C65-6188I-S                         | C65-6180I-S                          |