



Valves which will be installed in ovens or heated zones require a standoff assembly, which locates the actuator out of the heated zone and supports both the valve and the handle or actuator. The 5/8" outside diameter standoff tube extends through the oven wall and is secured by means of a clamp ring supplied with the assembly.

If you are converting an actuated valve from a closemount to a standoff application, order the appropriate clamp ring and two screws in addition to the standoff assembly. Consult the factory for availability of non-standard lengths.

The microelectric actuator for selectors uses a special standoff assembly (SOMMP) which is keyed to both valve and actuator. The key guarantees proper alignment and positioning of the valve.

Product numbers show the most common length of standoffs: 4" for air actuators and manual knobs, 2" for microelectric and standard electric actuators. Standoff assemblies are available in lengths of 2", 3", and 6". To order a 6" standoff instead of a 4" one, change the 4 in the product number to a 6.

### MORE INFORMATION

For illustrations of standoffs on valves and actuators, see pages 206-207.

### TECH TIP

If you need the **actuator as well as the hardware**, you can order it complete with the appropriate hardware or with the required standoff already installed.

### Actuators

Air ..... pages 194-195  
Microelectric .... 188-191  
Standard elec ..... 193

### Standoff assemblies and mounting hardware

*for actuators*

	Standoff assembly <i>Prod No</i>	Clamp ring <i>Prod No</i>	Screws <i>Prod No</i>
<b>Air actuators</b>			
For Valco two position valves with 1 or 2 mounting holes	450A	CR3	HWSC-SC8-6
with no mounting holes	450AMP	CR3	HWSC-SC8-6
For Valco selectors	450AMP	CR3	HWSC-SC8-6
For Cheminert valves	450AMP	CR3	HWSC-SC8-6
<b>Microelectric actuators</b>			
For Valco two position valves with 1 or 2 mounting holes	250A	CR8	HWSC-SC8-8B
with no mounting holes	250AMP	CR8	HWSC-SC8-8B
For Valco multiposition valves (UW and MW Types only)	250AMMP	CR10	HWSC-SC8-8TDH
For Cheminert two position valves	250AMP	CR8	HWSC-SC8-8B
For Cheminert selectors	250AMMP	CR10	HWSC-SC8-8TDH
<b>Standard electric actuators</b>			
For Valco two position valves with 1 or 2 mounting holes	250A	CR3	HWSC-SC8-8B
with no mounting holes	250AMP	CR3	HWSC-SC8-8B
For Valco selectors	250AMP	CR3	HWSC-SC8-8B
For Cheminert valves	250AMP	CR3	HWSC-SC8-8B

### Standoff assemblies

*for manual valves*

Includes knob, standoff assembly, retainer, and adapter. For illustration, see page 206, top.

	<i>Prod No</i>
For Valco W and UW Type two position valves rated less than 5,000 psi with 1 or 2 mounting holes	450WK
with no mounting holes	450WKMP
For Valco UW Type two position valves rated 5,000 psi and higher with 1 or 2 mounting holes	450UTH
with no mounting holes	450UTHMP
For Cheminert valves	450WKMP

