

On/off valves

SPECS			Manual		Manual with 4" standoff	Air actuated with 1" standoff	Air actuated with 4" standoff
Temp	Pressure		Fitting size	Bore	Prod No	Prod No	Prod No
Standard temperature			Standard temperature				
1/16"	100°C	10,000 psi	1/16"	0.50 mm	SFVO	-	ASFVO
				0.75 mm	SFVOL	-	ASFVOL
High temperature / high pressure			High temperature / high pressure				
1/16"	300°C	6,000 psi	1/16"	0.50 mm	SFVOHT	SFVOHT4	ASFVOHT
1/8"	300°C	2,000 psi		0.75 mm	-	-	ASFVOLHT
			1/8"	1.50 mm	-	-	ASFVO2HT

Prime/purge valves

SPECS			Manual		Air actuated with 1" standoff	Air actuated with 4" standoff
Temp	Pressure		Fitting size	Bore	Prod No	Prod No
Standard temperature			Standard temperature			
1/16"	100°C	10,000 psi	1/16"	0.50 mm	SFV	ASFV
				0.75 mm	SFVL	ASFVL
High temperature / high pressure			High temperature / high pressure			
1/16"	300°C	6,000 psi	1/16"	0.50 mm	-	ASFVHT
1/8"	300°C	2,000 psi		0.75 mm	-	ASFVLHT
			1/8"	1.50 mm	-	ASFV2HT

