

ValcoBond Columns

VB-1

100% dimethylpolysiloxane

REPLACES

DB-1, DB-1ms, HP-1, HP-1MS, Ultra-1, Rtx-1, Rtx-1MS, SPB-1, MDN-1, BP-1, CP-Sil 5 CB, GB-1, 007-1, OV-1, SE-30, AT-1 and ZB-1

0.10 mm ID	df	Prod No
10 meters	0.10	CFS-A01010-010B
10 meters	0.20	CFS-A01010-020B
10 meters	0.40	CFS-A01010-040B
20 meters	0.10	CFS-A02010-010B
20 meters	0.20	CFS-A02010-020B
20 meters	0.40	CFS-A02010-040B

0.18 mm ID	df	Prod No
10 meters	0.10	CFS-A01018-010B
10 meters	0.18	CFS-A01018-018B
10 meters	0.40	CFS-A01018-040B
30 meters	0.10	CFS-A03018-010B
30 meters	0.18	CFS-A03018-018B
30 meters	0.40	CFS-A03018-040B
30 meters	1.00	CFS-A03018-100B
40 meters	0.18	CFS-A04018-018B
40 meters	0.40	CFS-A04018-040B

0.25 mm ID	df	Prod No
15 meters	0.10	CFS-A01525-010B
15 meters	0.25	CFS-A01525-025B
15 meters	0.50	CFS-A01525-050B
15 meters	1.00	CFS-A01525-100B
15 meters	1.50	CFS-A01525-150B
30 meters	0.10	CFS-A03025-010B
30 meters	0.25	CFS-A03025-025B
30 meters	0.50	CFS-A03025-050B
30 meters	1.00	CFS-A03025-100B
30 meters	1.50	CFS-A03025-150B
60 meters	0.10	CFS-A06025-010B
60 meters	0.25	CFS-A06025-025B
60 meters	0.50	CFS-A06025-050B
60 meters	1.00	CFS-A06025-100B
60 meters	1.50	CFS-A06025-150B

0.32 mm ID	df	Prod No
15 meters	0.10	CFS-A01532-010B
15 meters	0.25	CFS-A01532-025B
15 meters	0.32	CFS-A01532-032B
15 meters	0.50	CFS-A01532-050B
15 meters	1.00	CFS-A01532-100B

0.32 mm ID, cont'd	df	Prod No
15 meters	2.00	CFS-A01532-200B
15 meters	3.00	CFS-A01532-300B
15 meters	4.00	CFS-A01532-400B
15 meters	5.00	CFS-A01532-500B
30 meters	0.10	CFS-A03032-010B
30 meters	0.25	CFS-A03032-025B
30 meters	0.32	CFS-A03032-032B
30 meters	0.50	CFS-A03032-050B
30 meters	1.00	CFS-A03032-100B
30 meters	2.00	CFS-A03032-200B
30 meters	3.00	CFS-A03032-300B
30 meters	4.00	CFS-A03032-400B
30 meters	5.00	CFS-A03032-500B
60 meters	0.10	CFS-A06032-010B
60 meters	0.25	CFS-A06032-025B
60 meters	0.32	CFS-A06032-032B
60 meters	0.50	CFS-A06032-050B
60 meters	1.00	CFS-A06032-100B
60 meters	2.00	CFS-A06032-200B
60 meters	3.00	CFS-A06032-300B
60 meters	4.00	CFS-A06032-400B
60 meters	5.00	CFS-A06032-500B

0.53 mm ID	df	Prod No
15 meters	0.15	CFS-A01553-015B
15 meters	0.50	CFS-A01553-050B
15 meters	1.00	CFS-A01553-100B
15 meters	1.50	CFS-A01553-150B
15 meters	3.00	CFS-A01553-300B
15 meters	5.00	CFS-A01553-500B
30 meters	0.15	CFS-A03053-015B
30 meters	0.50	CFS-A03053-050B
30 meters	1.00	CFS-A03053-100B
30 meters	1.50	CFS-A03053-150B
30 meters	3.00	CFS-A03053-300B
30 meters	5.00	CFS-A03053-500B
60 meters	1.00	CFS-A06053-100B
60 meters	1.50	CFS-A06053-150B
60 meters	3.00	CFS-A06053-300B
60 meters	5.00	CFS-A06053-500B

PRIMARY APPLICATIONS

Amines
Flavors
Fragrances
Hydrocarbons
Pesticides
PCBs
Phenols
Sulfur compounds
EPA Methods
504, 551, 1618
NIOSH Methods
1300-1301,
1400-1403,
1450, 1501, 2005

TO ORDER

For prices or more information about your specific application, contact VICI Metronics:

Toll-free 877-737-1887
Tel360-697-9199
Fax360-697-6682

columns@vici.com

VB-35

(35%Phenyl)-methylpolysiloxane

REPLACES

DB-35, AT-35, MDN-35, DB-35ms, Rtx-35, BP-35, HP-35, Rtx-35MS, 007-11, HP-35MS, Sup-Herb, ZB-35

0.25 mm ID	df	Prod No
15 meters	0.25	CFS-C01525-025B
15 meters	0.50	CFS-C01525-050B
30 meters	0.25	CFS-C03025-025B
30 meters	0.50	CFS-C03025-050B
60 meters	0.25	CFS-C06025-025B
60 meters	0.50	CFS-C06025-050B

0.32 mm ID	df	Prod No
15 meters	0.25	CFS-C01532-025B
15 meters	0.50	CFS-C01532-050B
30 meters	0.25	CFS-C03032-025B
30 meters	0.50	CFS-C03032-050B

0.32 mm ID, cont'd	df	Prod No
60 meters	0.25	CFS-C06032-025B
60 meters	0.50	CFS-C06032-050B

0.53 mm ID	df	Prod No
15 meters	0.50	CFS-C01553-050B
15 meters	1.00	CFS-C01553-100B
30 meters	0.50	CFS-C03053-050B
30 meters	1.00	CFS-C03053-100B
60 meters	0.50	CFS-C06053-050B
60 meters	1.00	CFS-C06053-100B

PRIMARY APPLICATIONS

Drugs
Pesticides
Herbicides
PAHs
Pharmaceuticals
PCBs
EPA Method 8081A
(organochlorine pesticides)