

VB-5

(5% Phenyl)-methylpolysiloxane

PRIMARY APPLICATIONS

Drugs
Herbicides
Hydrocarbons
PCBs
Pesticides
Phenols
Semi-volatiles
Sulfur compounds

0.10 mm ID	df	Prod No
10 meters	0.10	CFS-B01010-010B
10 meters	0.20	CFS-B01010-020B
20 meters	0.10	CFS-B02010-010B
20 meters	0.20	CFS-B02010-020B
0.18 mm ID	df	Prod No
10 meters	0.18	CFS-B01018-018B
10 meters	0.40	CFS-B01018-040B
20 meters	0.18	CFS-B02018-018B
20 meters	0.40	CFS-B02018-040B
30 meters	0.18	CFS-B03018-018B
30 meters	0.40	CFS-B03018-040B
40 meters	0.18	CFS-B04018-018B
40 meters	0.40	CFS-B04018-040B
0.25 mm ID	df	Prod No
15 meters	0.10	CFS-B01525-010B
15 meters	0.25	CFS-B01525-025B
15 meters	0.50	CFS-B01525-050B
15 meters	1.00	CFS-B01525-100B
30 meters	0.10	CFS-B03025-010B
30 meters	0.25	CFS-B03025-025B
30 meters	0.50	CFS-B03025-050B
30 meters	1.00	CFS-B03025-100B
60 meters	0.10	CFS-B06025-010B
60 meters	0.25	CFS-B06025-025B
60 meters	0.50	CFS-B06025-050B
60 meters	1.00	CFS-B06025-100B
0.32 mm ID	df	Prod No
15 meters	0.10	CFS-B01532-010B
15 meters	0.25	CFS-B01532-025B
15 meters	0.50	CFS-B01532-050B
15 meters	1.00	CFS-B01532-100B
15 meters	2.00	CFS-B01532-200B
15 meters	3.00	CFS-B01532-300B
15 meters	5.00	CFS-B01532-500B

0.32 mm ID, cont'd	Prod No	
30 meters	0.10	CFS-B03032-010B
30 meters	0.25	CFS-B03032-025B
30 meters	0.50	CFS-B03032-050B
30 meters	1.00	CFS-B03032-100B
30 meters	2.00	CFS-B03032-200B
30 meters	3.00	CFS-B03032-300B
30 meters	5.00	CFS-B03032-500B
0.53 mm ID	df	Prod No
15 meters	0.50	CFS-B01553-050B
15 meters	1.00	CFS-B01553-100B
15 meters	1.50	CFS-B01553-150B
15 meters	2.00	CFS-B01553-200B
15 meters	2.65	CFS-B01553-265B
15 meters	3.00	CFS-B01553-300B
15 meters	5.00	CFS-B01553-500B
30 meters	0.50	CFS-B03053-050B
30 meters	1.00	CFS-B03053-100B
30 meters	1.50	CFS-B03053-150B
30 meters	2.65	CFS-B03053-265B
30 meters	3.00	CFS-B03053-300B
30 meters	5.00	CFS-B03053-500B
60 meters	1.00	CFS-B06053-100B
60 meters	1.50	CFS-B06053-150B
60 meters	2.00	CFS-B06053-200B
60 meters	2.65	CFS-B06053-265B
60 meters	3.00	CFS-B06053-300B
60 meters	5.00	CFS-B06053-500B

REPLACES

DB-5, DB-5ms,
HP-5, HP-5MS,
Ultra-5, Rtx-5, Rtx-
5MS, Rtx-5sil MS,
SPB-5, MDN-5,
BP-5, CP-Sil 8 CB,
GB-5, 007-5, OV-5,
SE-54, AT-5, and
ZB-5

VB-50/608

(50%Phenyl)-methylpolysiloxane

PRIMARY APPLICATIONS

Drugs
Pharmaceuticals
Herbicides
Steroids
PAHs
Tocopherols
PCBs
EPA Methods
Pesticides
508, 608 and 8080

0.25 mm ID	df	Prod No
15 meters	0.25	CFS-D01525-025B
15 meters	0.50	CFS-D01525-050B
30 meters	0.15	CFS-D03025-015B
30 meters	0.25	CFS-D03025-025B
30 meters	0.50	CFS-D03025-050B
60 meters	0.25	CFS-D06025-025B
60 meters	0.50	CFS-D06025-050B
0.32 mm ID	df	Prod No
15 meters	0.25	CFS-D01532-025B
15 meters	0.50	CFS-D01532-050B
15 meters	1.00	CFS-D01532-100B
30 meters	0.25	CFS-D03032-025B
30 meters	0.50	CFS-D03032-050B
30 meters	1.00	CFS-D03032-100B

0.32 mm ID, cont'd	Prod No	
60 meters	0.25	CFS-D06032-025B
60 meters	0.50	CFS-D06032-050B
60 meters	1.00	CFS-D06032-100B
0.53 mm ID	df	Prod No
15 meters	0.50	CFS-D01553-050B
15 meters	0.83	CFS-D01553-083B
15 meters	1.00	CFS-D01553-100B
30 meters	0.50	CFS-D03053-050B
30 meters	0.83	CFS-D03053-083B
30 meters	1.00	CFS-D03053-100B
60 meters	0.50	CFS-D06053-050B
60 meters	0.83	CFS-D06053-083B
60 meters	1.00	CFS-D06053-100B

REPLACES

DB-17, AT-50,
SP-2250, DB-17ms,
BPX-50, SP-17,
DB-608, 007-17,
SPB-608, HP-50+,
SPB-50, ZB-50,
Rtx-50