

GC CAPILLARY COLUMNS from VICI Metronics

ValcoBond Columns

VB-Wax

100% bonded polyethylene glycol

REPLACES	0.10 mm ID	df	Prod No	0.32 mm ID	df	Prod No
DB-WAX, DB-WAXetr, HP-WAX, HP-InnoWAX, HP-20M, CB-WAX, Stabilwax, RtxWAX, SUPEROX II, SUPELCOWAX-10, BP-20, CP-WAX, 52 CB, GB-WAX, 007-CW, OV-WAX, AT-WAX, and ZB-WAX	10 meters	0.10	CFS-G01010-010A	15 meters	0.25	CFS-G01532-025A
	20 meters	0.10	CFS-G02010-010A	15 meters	0.50	CFS-G01532-050A
				15 meters	1.00	CFS-G01532-100A
				30 meters	0.25	CFS-G03032-025A
				30 meters	0.50	CFS-G03032-050A
				30 meters	1.00	CFS-G03032-100A
				60 meters	0.25	CFS-G06032-025A
				60 meters	0.50	CFS-G06032-050A

PRIMARY APPLICATIONS

Alcohols
Aldehydes
Aromatics
Flavors
Fragrances
Organic Acids
Solvents

VB-624/1301

(6% Cyanopropyl-phenyl)-methylpolysiloxane

REPLACES	0.18 mm ID	df	Prod No	0.32 mm ID	df	Prod No
DB-624, HP-624, HP-VOC, Rtx-624, Rtx-Volatiles, BP-624, Vocol, 007-624, 007-502, NON-PAKD, 624, ZB-624	10 meters	1.00	CFS-E01018-100A	15 meters	1.80	CFS-E01532-180A
	20 meters	1.00	CFS-E02018-100A	30 meters	1.80	CFS-E03032-180A
	30 meters	1.00	CFS-E03018-100A	60 meters	1.80	CFS-E06032-180A
	40 meters	1.00	CFS-E04018-100A			

PRIMARY APPLICATIONS

EPA Methods
501.3
502.2
503.1
524.2
601
602
8010
8015
8020
8240

VB-1701

(14% Cyanopropyl-phenyl)-methylpolysiloxane

REPLACES	0.25 mm ID	df	Prod No	0.32 mm ID, cont'd	Prod No	
DB-1701, 007-1701, HP-1701, CP-Sil 19 CB, Rtx-1701, SPB-1701, BP-10, ZB-1701	15 meters	0.25	CFS-F01525-025A	60 meters	0.25	CFS-F06032-025A
	15 meters	0.50	CFS-F01525-050A	60 meters	0.50	CFS-F06032-050A
	30 meters	0.25	CFS-F03025-025A	60 meters	1.00	CFS-F06032-100A
	30 meters	0.50	CFS-F03025-050A			

PRIMARY APPLICATIONS

Drugs
PAHs
PCBs
Pesticides
Phenols
Solvents
Tranquilizers

TO ORDER

Toll-free 877-737-1887
Tel360-697-9199
Fax.....360-697-6682

columns@vici.com

CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034