

ValcoPLOT Columns

Molesieve 5Å

Molesieve 5Å

REPLACES

GS-Molesieve 5A
HP-PLOT Molesieve
CP-Molesieve 5A
Rt-Msieve-5A
MXT-Msieve-51
PLT-5A

ValcoPLOT Molesieve 5Å PLOT columns offer greatly enhanced analytical efficiency at economical prices. Our thick film columns separate Ar/O₂ without the need for cryogenic equipment. ValcoPLOT Molesieve 5Å PLOT thin film columns offer fast elution of carbon monoxide with near perfect peak symmetry, and our proprietary bonding technology ensures that the particles stay put even when columns are used with valves.

PRIMARY APPLICATIONS
Gases

TO ORDER

For prices or more information about your specific application, contact VICI Metronics:

Toll-free 877-737-1887
Tel360-697-9199
Fax360-697-6682

columns@vici.com


Fused silica

0.53 mm ID	df	Prod No
15 meters	20	CFS-X1553-200
15 meters	50	CFS-X1553-500
30 meters	20	CFS-X3053-200
30 meters	50	CFS-X3053-500

Stainless steel

0.53 mm ID	df (µm)	Prod No
15 meters	20	CSS-X1553-200
30 meters	20	CSS-X3053-200
30 meters	50	CSS-X3053-500


ValcoPLOT Molesieve 5Å – Stainless steel



MOLESIEVE TEST MIX

Column: Stainless steel
Molesieve PLOT
30 m x 0.53 mm x 20.0 µm
Column temp: 40°C isothermal
Detector: PDD Model D-4-I

ValcoPLOT Molesieve 5Å – Fused silica



MOLESIEVE TEST MIX

Column: Fused silica
Molesieve PLOT
30 m x 0.53 mm x 20.0 µm
Column temp: 40°C isothermal
Detector: PDD Model D-4-I

CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034