



Analytical Syringes, Valves, Probes, and Custom Bent Tubing from VICI Precision Sampling

Micro Valves for GC and LC

- 200 psi, .060" bore
- Compact 1" design
- Convenient panel mount
- Variety of configurations

Simplify your liquid or gas handling application with a VICI Precision Sampling Micro valve. The unique design of the fitting detail allows a leak-free seal with no potential for rotor damage from overtightening. Internal parts are PEEK and PTFE.



(Fittings not included.
For example fittings, see page 68.)

Micro valves for GC and LC

1/4-28

	Prod No	SPECS		
"T" flow path	3 ports	660100	3 PORT	"T"
	4 ports	660110	4 PORT	
180° flow path	2 ports	660200	2 PORT	180°
	4 ports	660210	4 PORT	
90° flow path	2 ports	660300	2 PORT	90°
	3 ports	660310	3 PORT	
	4 ports	660320	4 PORT	

MORE INFORMATION
1/4-28 fittings . . . pp 68-78

TO ORDER
Toll-free 800-828-1653
Tel225-927-1128
Toll-free
fax866-429-7741
Fax225-923-1331

FOR OUR COMPLETE LINE OF PRODUCTS
Visit our website at
viciprecisionsampling.com
or call us for a catalog:

CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034