

## Decoding Valco Selector Product Numbers

Examples:

**A 2 VL SC 6 MW E2**

**A2VLSC6MWE2:** Air actuated, 2" standoff, 1/4" valve, SC flowpath, 6 positions, MW type, Valcon E2 rotor, standard Nitronic 60 body

**EMT 4 C SD 4 UW**

**EMT4CSD4UW:** Microelectric actuator, 4" standoff, 1/16" valve, SD flowpath, 4 positions, UW type, Valco H (blank = H) rotor, standard N60 body

**E 3 ST 10 MW T HC**

**E3ST10MW5HC:** Standard electric actuator, 3" standoff, 1/8" (blank = 1/8") valve, ST flowpath, 10 positions, MW type, Valcon T rotor, Hastelloy C body

NUMBER OF POSITIONS	VALVE TYPE	ROTOR MATERIAL	SPECIAL BODY MATERIAL
6. Required.	7. Required.	8. Required.	9. Optional. Body material is
6	Low pressure	E Valcon E	HC Hastelloy C
8	UW high pressure	E2 Valcon E2	IN Inconel 600
10		M Valcon M	M4 Monel 400
12		P Valcon P	NI Nickel 200
16		R Valcon R	N5 Nitronic 50
		T Valcon T	TI Titanium
		TF Valcon TF	

**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd  
Australian Distributors; Importers & Manufacturers

**NOTE!**

This chart is for decoding existing product numbers, **not** for inventing new ones. Some options can not work with certain valve types and designs!