



**New Features**

Calculates linear velocity based on column ID.  
Includes standard RS232 ComPort with printer function.

Has auto zero function,  
Includes rechargeable battery pack and charger.

**Also**

Reads flow accurately to 500ml/min.

Measures N2, He, H., 5% Ar/Me, and Air.

Measures split flow and mass flow.

Has pulse free operation that will not interfere with EPCs.



Includes auto "off" battery saver.

Includes NIST traceable calibration certificate.

Includes capillary column adaptor.

**Restek's Veri-Flow 500 Electronic Flow Meter**

110V, cat.# 21643 each  
220V, cat.# 21645 each



- Linear velocity measured based on column ID.
- Linear flow range from 0.05 to 20ml/min. (He).
- NIST traceable for He & H2.
- Manual zero function.
- Battery operated for easy portability.
- Capillary column adaptor.

**Restek's MicroFlo 20 Flow Meter**

cat.# 21644 each