



Gold Inlet Seals for HP GCs

Deactivated.
Prepacked with deactivated fused silica wool.
Purchased individually and in 5- or 25-packs.
Guaranteed to meet original equipment specifications.

Shimadzu 17A PTV Sleeve

cat.# 21705, each
cat.# 21706, 5-pack
cat.# 21707, 25-pack

**COMING SOON~ Inlet Sleeves for Varian 1078
instruments. Call your local distributor for details!**

Red or black discoloration has occurred sporadically with some Restek Gold Inlet Seals. Initially this problem was hard to diagnose and was believed to be amplified by certain sample contaminants. Restek's chemists investigated this phenomena and came up with a new plating process to prevent discoloration. By using this new plating technology, we have eliminated any discoloration. Please call Restek's technical service at 800-356-1688, ext. 4, if you have questions about these improved inlet seals.

0.8mm ID Opening:
cat.# 21317,2-pk; cat.# 21318, 10-pk.
1.2mm ID Opening:
cat.# 21305,2-pk; cat.# 21306, 10-pk.

Restek feels a strong responsibility toward helping analysts keep abreast of all the latest advancements by providing free literature. Our literature provides late-breaking information on capillary technology.

Free Technical Literature / Catalogs:

- 1996 Chromatography Products
- Chemical Standards
- HPLC Columns and Accessories
- Packed Columns
- MXT Columns and Accessories
- Clinical, Forensic, & Pharmaceutical Chromatography Products
- Performance Products for HP GCs
- Air Monitoring Products

Technical Guides:

Capillary Column Installation Guide
A Guide to GC Setup

- Hints for Analyzing Volatile Organics
- A Guide to Mirzimirg Septa Problems
- Operating Hints for Split/Splitless Injection
- A Guide When Injecting Dir8 Samples
- Guide to Direct/On-Column Flash Vaporization Injection
- Guide to the Analysis of Petroleum Hydrocarbons in Soil & Water
- Capillaq Column Installation Guide

Other Literature Available:

- The Restek Advantage newsletter
- Application Notes and Technical Tips (covering environmental, petrochemical, clinical/forensic, and food/flavor analyses, as well as new capillary columns and chromatography techniques)
- Capillary Chromatography Essentials wall chart
- Sample Extraction Products
- Silcosteel® Products
- SilcoCan™ Air Sampling Canisters
- “Environmental Pinnacle” HPLC Columns

Restek (US)

110 Benner Circle
Bellefonte, PA 16823
Phone: (800) 356-1688
FAX: (814) 353-1309

Restek GmbH

Sulzbacher Str. 15
D-658 12 Bad Soden, Germany
Telefon: 06196/65130
Telefax: 06196/623 01

Restek France

1, Rue Montespan
91024 Evry, France
Telefon: (1) 60.78.32.10
Telefax: (1) 60.78.79.90

Restek Trademarks: CarboFrit, Integra-Guard, Resprep, Diskcover, Rtx, MXT Silcostee I, SilcoCan, Pinnacle, RedLite, Thermolite, and the Restek logo.

Other Trademarks: Empore (3M Company), Thermogreen (Supelco, Inc.), Chromseal 9001 HT (Chromacol), Teflon (E.I. du Pont de Nemours inc.).



MasterCard®, VISA®, and American Express® are accepted.