



Sturdy Teflon" valve offers inertness and longevig.

**RESPREP-12T Extraction System for SPE Cartridges: cat.# 24001**

Complete Kit includes: propylene top with twelve sturdy Teflona flow regulation valves, 1/16" Teflon" tubing, glass block with built in vacuum regulator and scratch resistant polypropylene base, multipurpose sample holding rack, and convenient valve turning tool.

*An Application Note is available detailing the Florisil" cleanup procedure for Organochlorine Pesticides and PCBs via the USEPA CLP methodology. Call your local distributor to receive a copy.*

Name	Sorbent Mass (mg)	Cartridg Volume(ml)	Frit Style*	Pkg. Size	Catalog Number
Florisil" Cartridges	500	3	PE	50	24031
Florisil" Cartridges	500	3	s s	50	24032
Florisil" Cartridges	500	10	PE	50	24033
Florisil" Cartridges	1000	6	PE	30	24034
Silica Gel Cartridges	500	3	PE	50	24035
Silica Gel Cartridges	500	3	s s	50	24034
Silica Gel Cartridges	500	10	PE	50	24037
Silica Gel Cartridges	1000	6	PE	30	24038
Alumina-N Cartridges	500	3	PE	50	24039
Alumina-N Cartridges	500	3	s s	50	24040
Alumina-N Cartridges	500	10	PE	50	24041
Alumina-A Cartridges	500	3	PE	50	24042
Alumina-A Cartridges	500	3	s s	50	24043
Alumina-A Cartridges	500	10	PE	50	24044
Alumina-B Cartridges	500	3	PE	50	24045
Alumina-B Cartridges	500	3	s s	50	24046
Alumina-B Cartridges	500	10	PE	50	24047

\*PE: Polyethylene frits; X3: Stain&s steel frits