


Peak Performers Silcosteel@ Fittings & Tubing


We've reduced Silcosteel@ fitting costs and are passing the savings on to you!

Don't see what you need? Send us your parts and we'll treat them for you! Remember, call first to obtain an authorization number!



Unions


sizes	cat.#
1/16"	20510
1/8"	20511
1/4"	20512



Elbows


sizes	cat.#
1/16"	20516
1/8"	20517
1/4"	20518

Each Silcosteel@ fitting (available in 1/16", 1/8", and 1/4-inch) is deactivated and lined with a fused silica layer to ensure a completely inert sample pathway.



Tees

sizes	cat.#
1/16"	20513
1/8"	20514
1/4"	20515



Reducing Fittings

sizes	cat.#
1/16to1/8"	20519
1/16to1/4"	20520
1/8to1/4	20521

FITTINGS ARE NOW

20%

LESS THAN ORIGINAL PRICE!

Silcosteel@ Fused Silica Lined Stainless Steel Tubing (Coiled)

	0.011" ID (0.28mm) 0.022" OD (0.56mm)	0.021" ID (0.53mm) 0.029" OD (0.74mm)	0.010" ID (0.25mm) 1/16" OD (1.59mm)
6ft.	20570	20563	20550
25ft.	20571	20564	20551
50ft.	20572	20565	20552
200ft.	20573	20566	20553
>400ft.	20574	20567	20554

	0.020" ID (0.51mm) 1/16" OD (1.59mm)	0.030" ID (0.76mm) 1/16" OD (1.59mm)	0.040" ID (1.02mm) 1/16" OD (1.59mm)
6ft.	20524	20530	20538
25ft.	20525	20531	20539
50ft.	20526	20532	20540
200ft.	20527	20533	20541
>400ft.	20528	20534	20542

	0.085" ID (2.16mm) 1/8" OD (3.18mm)	0.210" ID (5.33mm) 1/4" OD (6.35mm)
6ft.	20545	20555
25ft.	20546	20556
5ft.	20547	20557
200ft.	20548	20558
>400ft.	20549	20559

Call for availability of lengths greater than 400ft.

*metric conversion: 6ft. (1.8m), 25ft. (7.6m), 50ft. (15.2m), 200ft. (61m), >400ft. (> 122m)

Our Silcosteel@ products are used by:

Hewlett-Packard, Varian, NASA, DOE, Shimadzu, Tekmar, OI, Entech, Graseby, Thermedics, Applied Automation, and ABBProcess Analytics.

Isn't it time you tried SilcosteelP?

Restek Corporation