

Peak Performers



new! Inlet Liners for Varian 1078 Injectors

- Deactivated.
- Available for split or splitless analysis.
- Offered individually or in economical 5-packs.
- Guaranteed to meet original equipment specifications.

SPLIT LINER (similar to Varian part # 03-918464-00).

cat.# 21708, each

cat.# 21709, 5-pack



SPLITLESS LINER (similar to Varian part # 03-918466-00).

cat.# 21711, each

cat.# 21712, 5-pack




new! Capillary Grade Hydrocarbon Trap

- Removes trace carrier gas impurities.
- Purged with ultra high purity helium.
- Packed with high surface area activated carbon.
- Reduces baseline stabilization time after installation.

- 1/8" CAPILLARY GRADE HYDROCARBON TRAP:
cat.# 21991, each

- 1/4" CAPILLARY GRADE HYDROCARBON TRAP:
cat.# 21992, each



new! High Capacity Indicating Oxygen Trap

- Indicates color change as O₂ and H₂O are trapped.
- Reduces O₂ to less than 1ppb.
- Usable with all carrier gases.
- Three times the capacity of low capacity indicating traps.
- Can be regenerated.

1/8" HIGH CAPACITY INDICATING OXYGEN TRAP:
cat.# 21995, each


1/4" HIGH CAPACITY INDICATING OXYGEN TRAP:
cat.# 21996, each

Return unbroken Restek traps and give them a new life. Please call our customer service group at 800-356-1688, ext. 4, for an RA number before returning.

REGENERATED TRAP:

1/8" fitting: cat.# 21993, each

1/4" fitting: cat.# 21994, each



Restek (U.S.):
110 Benner Circle
Bellefonte, PA 16823
Phone:(800)356-688
FAX: (814) 353-1309
<http://www.restekcorp.com>

Restek GmbH:
Sulzbacher Str. 15
D-658 12 Bad Soden, Germany
Telefon: 06196165130
Telefax: 06 19616230 1

Restek Trademarks: Rtx, Silcosteel, SilcoCan, SilcoPort, Stabilwax, Crossbond, Press;
Other trademarks: HP-5MS, DB-5MS, SE-30, OV-101, Chromosorb, Carbowax, Ultra-Seal