



Space Saving, Personal-Sized Centrifuge!

**Compact mini-centrifuge is ideal** for microfiltration and quick spin-downs of microcentrifuge tubes.

**Patented Auto Lid opens and closes automatically** with one-finger touch.

**Easy-to-use dual start/stop function** – start and stop using power button or by opening and closing the lid.

**Rubber footing keeps unit in place** and minimizes vibration.

**Easy to clean** ABS body and rotor and clear polycarbonate cover.

**Removable cord** for easy storage.



**Complete set includes** rotors for 1.5ml and 2.0ml tubes, adapters for .5ml tubes, and a strip rotor to accommodate two 8 x .2ml tubes.



Ordering Information

Catalog No.	Description	Price
		<i>Special</i>
PC-200	Whale Mini-Centrifuge	A\$400.00 AUD=0.80USD

Specifications

Supplied with	1.5/2.0ml rotor, .5ml adapters, and two 8-strip 0.2ml tube rotor	
Maximum Speed	5000 RPM/2000 x g	
Maximum Capacity		
Standard Tube Rotor	6 x 1.5/2.0ml tubes	 
Adapters	6 x .5ml tubes	
Strip Tube Rotor	2 x 8 x .2ml tubes	

Dimensions

Dimensions (D x W x H)	126mm x 138mm x 100mm (5.0 x 5.4 x 4.0in)
Weight	1.0 kg (2.2 lb)
Electrical	100-240V, 50/60 Hz (supplied with 120V cord)