

Miniature Air Sampling Canisters

- Ideal for indoor air, personal, emergency response, or soil gas sampling.
- 400cc or 1000cc.
- Low pressure applications not exceeding 40psig.
- Available with quick-connect fitting that is compatible with sampling and analysis instruments.
- Also available with non-treated or Sulfinert®-treated valve.

These small canisters are designed for controlled sampling, such as personal air sampling, as an alternative to tube and pump samplers. The 1000cc canister is suitable for sampling volatile organic compounds in air according to US EPA Methods TO-14 and TO-15.

Restek offers these products in stainless steel or with Sulfinert® treatment, for greatest inertness. We continue to offer passive coating technologies that are unmatched in the air sampling industry—try a Sulfinert®-treated canister and achieve the ultimate in analyte stability.

Miniature Air Sampling Canisters with Quick-Connect Stem Fittings

Description	Volume	qty.	cat.#
Electro-Polished Miniature Canister with Quick-Connect Stem Fitting	400cc	ea.	24188
	1000cc	ea.	24194
Sulfinert®-Treated Miniature Canister with Quick-Connect Stem Fitting	400cc	ea.	24189
	1000cc	ea.	24195
Sulfinert®-Treated Miniature Canister with Sulfinert®-Treated Quick-Connect Stem Fitting	400cc	ea.	24190
	1000cc	ea.	24196

Miniature Air Sampling Canisters with Metal-Seated Diaphragm Valve

Description	Volume	qty.	cat.#
Electro-Polished Miniature Canister with Metal-Seated Diaphragm Valve	400cc	ea.	24191
	1000cc	ea.	24197
Sulfinert®-Treated Miniature Canister with Metal-Seated Diaphragm Valve	400cc	ea.	24192
	1000cc	ea.	24198
Sulfinert®-Treated Miniature Canister with Sulfinert®-Treated Diaphragm Valve	400cc	ea.	24193
	1000cc	ea.	24199


Miniature Air Sampling Canisters with Nut & Ferrule **new!**

Description	Volume	qty.	cat.#
Electro-Polished Miniature Canister with Nut & Ferrule	400cc	ea.	24205
	1000cc	ea.	24206
Sulfinert®-Treated Miniature Canister with Nut & Ferrule	400cc	ea.	24207
	1000cc	ea.	24208

Quick-Connect Fittings for Miniature Air Sampling Canisters

Connection: 1/4" tube fitting.

Description	qty.	cat.#
Quick-Connect Stem Fitting	ea.	24185
Sulfinert®-Treated Quick-Connect Stem Fitting	ea.	24186
Quick-Connect Body Fitting	ea.	24187



Dimensions:
400cc = 2.75" diameter,
5.35" long
1000cc = 2.75" diameter,
11.92" long



a plus 1 story

"Restek has consistently provided high quality chromatography columns and supplies to me for well over a decade. Over the last two years, I have extensively been involved with air analysis, TO-15, etc., and Restek provides the highest quality canisters, mini-cans, and flow controllers in the market today."

Scott Van Etten, IH Laboratory Manager, EMSL Analytical

for
**reference
books**

see pages 450-459