


Sunniest Guard Cartridge



Photo of Cartridge and Holder




Schematic Diagram of Cartridge 10x4mm

Feature


- *Simple structure
- *Low dead volume

Comparison of chromatograms

Analytical column only



Guard and Analytical column



Column: Sunniest C18, 5 μ m 150 x 4.6 mm

Guard cartridge 10 x 4 mm

Mobile phase:

CH₃OH/20mM Phosphate buffer pH7.5 = 80/20

Flow rate: 1.0 mL/min

Temperature: 40 °C

Sample: 1 = Uracil,

2 = Propranolol,

3 = Nortriptyline,

4 = Amitriptyline,

TF: USP tailing factor

Ordering information

Guard cartridge (4-pak) 10 x 4 mm	Particle size	Price
Guard cartridge (4-pak) Aqua	5 μ m	JPY 25,200
Guard cartridge (4-pak) C18	5 μ m	JPY22,000
Guard cartridge holder 10 x 4 mm	---	JPY 12,000