







FTLL250	NUMBER	MATERIAL
Female Luer Thread Style to 200 Series Barb, 3/16" (4.8 mm) ID Tubing 	FTLL250-1	White Nylon
	FTLL250-3	Red Nylon
	FTLL250-5	Blue Nylon
	FTLL250-6	Natural Polypropylene
	FTLL250-9	Clear Polycarbonate
	FTLL250-9002	Rad. Stable Polycarbonate
	FTLL250-J1A	Natural Kynar PVDF


FTLL450	NUMBER	MATERIAL
Female Luer Thread Style to 400 Series Barb, 3/16" (4.8 mm) ID Tubing 	FTLL450-1	White Nylon
	FTLL450-6	Natural Polypropylene
	FTLL450-9	Clear Polycarbonate
	FTLL450-9002	Rad. Stable Polycarbonate
	FTLL450-J1A	Natural Kynar PVDF

FTLL055	NUMBER	MATERIAL
Female Luer Thread Style to 500 Series Barb, 1/4" (6.4 mm) ID Tubing 	FTLL055-1	White Nylon
	FTLL055-2	Black Nylon
	FTLL055-3	Red Nylon
	FTLL055-4	Green Nylon
	FTLL055-5	Blue Nylon
	FTLL055-06	Orange Nylon
	FTLL055-07	Yellow Nylon
	FTLL055-6	Natural Polypropylene
	FTLL055-9	Clear Polycarbonate
	FTLL055-9002	Rad. Stable Polycarbonate

FTLL360	NUMBER	MATERIAL
Female Luer Thread Style to 300 Series Barb, 1/4" (6.4 mm) ID Tubing 	FTLL360-1	White Nylon
	FTLL360-3	Red Nylon
	FTLL360-5	Blue Nylon
	FTLL360-6	Natural Polypropylene
	FTLL360-6005	Animal-Free Polypropylene
	FTLL360-9	Clear Polycarbonate
	FTLL360-9002	Rad. Stable Polycarbonate
FTLL360-J1A	Natural Kynar PVDF	

FTLLSB004	NUMBER	MATERIAL
Female Luer Thread Style Smooth Bore to 500 Series Barb, 1/16" (1.6 mm) ID Tubing 	FTLLSB004-1	White Nylon
	FTLLSB004-6	Natural Polypropylene
	FTLLSB004-9	Clear Polycarbonate

FTLLSB007	NUMBER	MATERIAL
Female Luer Thread Style Smooth Bore to 500 Series Barb, 3/32" (2.4 mm) ID Tubing 	FTLLSB007-1	White Nylon
	FTLLSB007-6	Natural Polypropylene
	FTLLSB007-9	Clear Polycarbonate

FTLLSB013	NUMBER	MATERIAL
Female Luer Thread Style Smooth Bore to 500 Series Barb, 1/8" (3.2 mm) ID Tubing 	FTLLSB013-1	White Nylon
	FTLLSB013-6	Natural Polypropylene
	FTLLSB013-9	Clear Polycarbonate

WHY CHOOSE A SMOOTH BORE?

Our Smooth Bore Luers have a smooth inner transition providing for easy insertion of guide wires, probes and other devices.

These components feature a specially designed finger grip to differentiate them from our other luer connectors.

TECHNICAL ASSISTANCE

Trying to choose the correct part? Need more material information than is published in this catalog? Call Value Plastics Technical Support for assistance at 970-267-5200.

ORDERING LARGE QUANTITIES?

Value Plastics offers quantity discounts on many of our fittings for orders over 5,000 pieces. Call 970-267-5200 for details

HRMalytic www.chromtech.net.au sales@chromtech.net.au
 ABN 14 643 445 058 **BY LID** Tel: (03) 9762 2034 Fax: +61 3 9761 1169
TECH nology ... for ALL Your chromatography supplies!
 Australian Distributors