

<b>MTLL430</b>	<b>NUMBER</b>	<b>MATERIAL</b>
Male Luer Integral Lock Ring to 400 Series Barb, 1/8" (3.2 mm) ID Tubing	MTLL430-1	White Nylon
	<b>MTLL430-3</b>	<b>Red Nylon</b>
	MTLL430-6	Natural Polypropylene
	MTLL430-9	Clear Polycarbonate
	MTLL430-9002	Rad. Stable Polycarbonate
	MTLL430-J1A	Natural Kynar PVDF



<b>MTLL250</b>	<b>NUMBER</b>	<b>MATERIAL</b>
Male Luer Integral Lock Ring to 200 Series Barb, 3/16" (4.8 mm) ID Tubing	MTLL250-1	White Nylon
	MTLL250-6	Natural Polypropylene
	<b>MTLL250-8003</b>	<b>Natural ABS</b>
	MTLL250-9	Clear Polycarbonate
	MTLL250-9002	Rad. Stable Polycarbonate
	MTLL250-J1A	Natural Kynar PVDF



<b>MTLL025</b>	<b>NUMBER</b>	<b>MATERIAL</b>
Male Luer Integral Lock Ring to 500 Series Barb, 5/32" (4.0 mm) ID Tubing	MTLL025-1	White Nylon
	MTLL025-2	Black Nylon
	MTLL025-3	Red Nylon
	<b>MTLL025-4</b>	<b>Green Nylon</b>
	MTLL025-5	Blue Nylon
	<b>MTLL025-06</b>	<b>Orange Nylon</b>
	<b>MTLL025-07</b>	<b>Yellow Nylon</b>
	MTLL025-6	Natural Polypropylene
	MTLL025-9	Clear Polycarbonate
	MTLL025-9002	Rad. Stable Polycarbonate



<b>MTLL450</b>	<b>NUMBER</b>	<b>MATERIAL</b>
Male Luer Integral Lock Ring to 400 Series Barb, 3/16" (4.8 mm) ID Tubing	MTLL450-1	White Nylon
	MTLL450-6	Natural Polypropylene
	MTLL450-9	Clear Polycarbonate
	<b>MTLL450-9002</b>	<b>Rad. Stable Polycarbonate</b>
	<b>MTLL450-J1A</b>	<b>Natural Kynar PVDF</b>



<b>MTLL240</b>	<b>NUMBER</b>	<b>MATERIAL</b>
Male Luer Integral Lock Ring to 200 Series Barb, 5/32" (4.0 mm) ID Tubing	MTLL240-1	White Nylon
	MTLL240-2	Black Nylon
	MTLL240-6	Natural Polypropylene
	MTLL240-9	Clear Polycarbonate
	MTLL240-9002	Rad. Stable Polycarbonate
	MTLL240-J1A	Natural Kynar PVDF



<b>MTLL445</b>	<b>NUMBER</b>	<b>MATERIAL</b>
Male Luer Integral Lock Ring to 400 Series Barb, .170" (4.3 mm) ID Tubing	MTLL445-1	White Nylon
	MTLL445-6	Natural Polypropylene
	MTLL445-9	Clear Polycarbonate
	<b>MTLL445-9002</b>	<b>Rad. Stable Polycarbonate</b>
	<b>MTLL445-J1A</b>	<b>Natural Kynar PVDF</b>



<b>MTLL055</b>	<b>NUMBER</b>	<b>MATERIAL</b>
Male Luer Integral Lock Ring to 500 Series Barb, 1/4" (6.4 mm) ID Tubing	MTLL055-1	White Nylon
	MTLL055-2	Black Nylon
	<b>MTLL055-3</b>	<b>Red Nylon</b>
	<b>MTLL055-4</b>	<b>Green Nylon</b>
	MTLL055-06	Orange Nylon
	MTLL055-07	Yellow Nylon
	MTLL055-6	Natural Polypropylene
	MTLL055-9	Clear Polycarbonate
	MTLL055-9002	Rad. Stable Polycarbonate



<b>MTLL035</b>	<b>NUMBER</b>	<b>MATERIAL</b>
Male Luer Integral Lock Ring to 500 Series Barb, 3/16" (4.8 mm) ID Tubing	MTLL035-1	White Nylon
	<b>MTLL035-2</b>	<b>Black Nylon</b>
	MTLL035-3	Red Nylon
	MTLL035-4	Green Nylon
	MTLL035-5	Blue Nylon
	<b>MTLL035-06</b>	<b>Orange Nylon</b>
	<b>MTLL035-07</b>	<b>Yellow Nylon</b>
	MTLL035-6	Natural Polypropylene
	MTLL035-9	Clear Polycarbonate
	MTLL035-9002	Rad. Stable Polycarbonate



<b>MTLL360</b>	<b>NUMBER</b>	<b>MATERIAL</b>
Male Luer Integral Lock Ring to 300 Series Barb, 1/4" (6.4 mm) ID Tubing	MTLL360-1	White Nylon
	MTLL360-6	Natural Polypropylene
	MTLL360-6005	Animal-Free Polypropylene
	MTLL360-9	Clear Polycarbonate
	MTLL360-9002	Rad. Stable Polycarbonate
	MTLL360-J1A	Natural Kynar PVDF

