


18210	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to 200 Series Barb, 1/16" (1.6 mm) ID Tubing	18210-1	White Nylon
	18210-2	Black Nylon
	18210-6	Natural Polypropylene
	18210-J1A	Natural Kynar PVDF



18410	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to 400 Series Barb, 1/16" (1.6 mm) ID Tubing	18410-1	White Nylon
	18410-2	Black Nylon
	18410-6	Natural Polypropylene
	18410-J1A	Natural Kynar PVDF



1810	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to Classic Series Barb, 1/16" (1.6 mm) ID Tubing	1810-1	White Nylon
	1810-2	Black Nylon
	1810-6	Natural Polypropylene
	1810-65	Powder Blue Polypropylene
	1810-J1A	Natural Kynar PVDF



18220	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to 200 Series Barb, 3/32" (2.4 mm) ID Tubing	18220-1	White Nylon
	18220-2	Black Nylon
	18220-6	Natural Polypropylene
	18220-J1A	Natural Kynar PVDF



18420	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to 400 Series Barb, 3/32" (2.4 mm) ID Tubing	18420-1	White Nylon
	18420-2	Black Nylon
	18420-6	Natural Polypropylene
	18420-J1A	Natural Kynar PVDF



1820	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to Classic Series Barb, 3/32" (2.4 mm) ID Tubing	1820-1	White Nylon
	1820-2	Black Nylon
	1820-6	Natural Polypropylene
	1820-63	Red Polypropylene
	1820-J1A	Natural Kynar PVDF



18230	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to 200 Series Barb, 1/8" (3.2 mm) ID Tubing	18230-1	White Nylon
	18230-2	Black Nylon
	18230-6	Natural Polypropylene
	18230-J1A	Natural Kynar PVDF



18430	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to 400 Series Barb, 1/8" (3.2 mm) ID Tubing	18430-1	White Nylon
	18430-2	Black Nylon
	18430-6	Natural Polypropylene
	18430-J1A	Natural Kynar PVDF



1830	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to Classic Series Barb, 1/8" (3.2 mm) ID Tubing	1830-1	White Nylon
	1830-2	Black Nylon
	1830-6	Natural Polypropylene
	1830-J1A	Natural Kynar PVDF



18240	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to 200 Series Barb, 5/32" (4.0 mm) ID Tubing	18240-1	White Nylon
	18240-2	Black Nylon
	18240-6	Natural Polypropylene
	18240-J1A	Natural Kynar PVDF



1840	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to Classic Series Barb, 5/32" (4.0 mm) ID Tubing	1840-1	White Nylon
	1840-2	Black Nylon
	1840-6	Natural Polypropylene
	1840-J1A	Natural Kynar PVDF



18445	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to 400 Series Barb, .170" (4.3 mm) ID Tubing	18445-1	White Nylon
	18445-2	Black Nylon
	18445-6	Natural Polypropylene
	18445-J1A	Natural Kynar PVDF



1850	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to Classic Series Barb, 3/16" (4.8 mm) ID Tubing	1850-1	White Nylon
	1850-2	Black Nylon
	1850-6	Natural Polypropylene
	1850-J1A	Natural Kynar PVDF



1860	NUMBER	MATERIAL
1/8-27 NPT Thread with 7/16" Hex to Classic Series Barb, 1/4" (6.4 mm) ID Tubing	1860-1	White Nylon
	1860-2	Black Nylon
	1860-6	Natural Polypropylene
	1860-J1A	Natural Kynar PVDF



HRomalytic www.chromotech.net.au sales@chromotech.net.au
 Australian Distributors: **TECH** www.chromotech.net.au
 ABN 14 643 445 058 **PT LTD** Tel: (03) 9762 2024 Fax: +61 3 9761 1169
 ... for ALL Your chromatography supplies !