

Material Suffix Key

-1 White Nylon	-10 Acetal	-2 Black Nylon	-3 Red Nylon
-4 Green Nylon	-5 Blue Nylon	-06 Orange Nylon	-07 Yellow Nylon
-6 Natural Polypropylene	-6005 Animal-Free Polypropylene	-8003 Natural ABS	-81 Lace White ABS
-9 Clear Polycarbonate	-9002 Radiation Stable Polycarbonate	-J1A Natural Kynar PVDF	

FTLL004 (-1 -2 -3 -4 -5 -06 -07 -6 -9 -9002 -J1A)

Female Luer Thread Style to 500 Series Barb, 1/16" (1.6 mm) ID Tubing


FTLL210 (-1 -2 -3 -5 -6 -9 -9002 -J1A)

Female Luer Thread Style to 200 Series Barb, 1/16" (1.6 mm) ID Tubing


FTLL410 (-1 -6 -9 -9002 -J1A)

Female Luer Thread Style to 400 Series Barb, 1/16" (1.6 mm) ID Tubing


FTLL007 (-1 -2 -3 -4 -5 -06 -07 -6 -9 -9002)

Female Luer Thread Style to 500 Series Barb, 3/32" (2.4 mm) ID Tubing


FTLL220 (-1 -3 -5 -6 -9 -9002 -J1A)

Female Luer Thread Style to 200 Series Barb, 3/32" (2.4 mm) ID Tubing


FTLL420 (-1 -6 -6005 -9 -9002 -J1A)

Female Luer Thread Style to 400 Series Barb, 3/32" (2.4 mm) ID Tubing


FTLL013 (-1 -2 -3 -4 -5 -06 -07 -6 -9 -9002)

Female Luer Thread Style to 500 Series Barb, 1/8" (3.2 mm) ID Tubing


FTLL230 (-1 -2 -3 -5 -6 -6005 -9 -9002 -J1A)

Female Luer Thread Style to 200 Series Barb, 1/8" (3.2 mm) ID Tubing


FTLL430 (-1 -6 -9 -9002 -J1A)

Female Luer Thread Style to 400 Series Barb, 1/8" (3.2 mm) ID Tubing


FTLL025 (-1 -2 -3 -4 -5 -06 -07 -6 -6005 -9 -9002)

Female Luer Thread Style to 500 Series Barb, 5/32" (4.0 mm) ID Tubing


FTLL240 (-0USP* -1 -2 -3 -5 -6 -9 -9002 -J1A)

Female Luer Thread Style to 200 Series Barb, 5/32" (4.0 mm) ID Tubing

*Natural Nylon USP Class VI


FTLL445 (-1 -6 -9 -9002 -J1A)

Female Luer Thread Style to 400 Series Barb, .170" (4.3 mm) ID Tubing


FTLL035 (-1 -2 -3 -4 -5 -06 -07 -6 -6005 -9 -9002 -J1A)

Female Luer Thread Style to 500 Series Barb, 3/16" (4.8 mm) ID Tubing


FTLL250 (-1 -3 -5 -6 -9 -9002 -J1A)

Female Luer Thread Style to 200 Series Barb, 3/16" (4.8 mm) ID Tubing


FTLL450 (-1 -6 -9 -9002 -J1A)

Female Luer Thread Style to 400 Series Barb, 3/16" (4.8 mm) ID Tubing


FTLL055 (-1 -2 -3 -4 -5 -06 -07 -6 -6005 -9 -9002)

Female Luer Thread Style to 500 Series Barb, 1/4" (6.4 mm) ID Tubing


FTLL360 (-1 -3 -5 -6 -9 -9002 -J1A)

Female Luer Thread Style to 300 Series Barb, 1/4" (6.4 mm) ID Tubing




FEMALE LUER:
Smooth Bore, Thread Style

FEMALE LUER:
Lug Style


Material Suffix Key

- | | | | |
|--------------------------|--------------------------------------|-------------------------|--------------------|
| -1 White Nylon | -10 Acetal | -2 Black Nylon | -3 Red Nylon |
| -4 Green Nylon | -5 Blue Nylon | -06 Orange Nylon | -07 Yellow Nylon |
| -6 Natural Polypropylene | -6005 Animal-Free Polypropylene | -8003 Natural ABS | -81 Lace White ABS |
| -9 Clear Polycarbonate | -9002 Radiation Stable Polycarbonate | -J1A Natural Kynar PVDF | |


FTLLSB004 (-1 -6 -9)
Female Luer Thread Style Smooth Bore to 500 Series Barb, 1/16" (1.6 mm) ID Tubing



FTLLSB007 (-1 -6 -9)
Female Luer Thread Style Smooth Bore to 500 Series Barb, 3/32" (2.4 mm) ID Tubing



FTLLSB013 (-1 -6 -9)
Female Luer Thread Style Smooth Bore to 500 Series Barb, 1/8" (3.2 mm) ID Tubing



WHY CHOOSE A SMOOTH BORE?

Our Smooth Bore Luers have a smooth inner transition providing for easy insertion of guide wires, probes and other devices.

These components feature a specially designed finger grip to differentiate them from our other luer connectors.

FEMALE LUER:
Lug Style

FTL10 (-1 -6 -9 -J1A)
Female Luer Lug Style to Classic Series Barb, 1/16" (1.6 mm) ID Tubing



FTL210 (-1 -6 -9 -J1A)
Female Luer Lug Style to 200 Series Barb, 1/16" (1.6 mm) ID Tubing



FTL20 (-1 -07 -6 -9 -J1A)
Female Luer Lug Style to Classic Series Barb, 3/32" (2.4 mm) ID Tubing



FTL220 (-1 -3 -5 -6 -9 -J1A)
Female Luer Lug Style to 200 Series Barb, 3/32" (2.4 mm) ID Tubing



FTL30 (-1 -6 -9 -J1A)
Female Luer Lug Style to Classic Series Barb, 1/8" (3.2 mm) ID Tubing



FTL230 (-1 -3 -5 -6 -9 -J1A)
Female Luer Lug Style to 200 Series Barb, 1/8" (3.2 mm) ID Tubing



FTL240 (-1 -6 -9)
Female Luer Lug Style to 200 Series Barb, 5/32" (4.0 mm) ID Tubing



LUER:
Adapters, Couplers and Lock Rings

FSLLR (-1 -2 -3 -4 -5 -06 -07 -6 -9 -9002)
Finger Snap Luer Lock Ring



MLRLC (-1 -6 -6005 -9)
Male Luer Coupler (May be used with separate rotating lock rings; FSLLR)



MTLCS (-1 -6 -9 -J1A)
Male Luer Slip Coupler



FTLLC (-1 -6 -9 -9002 -J1A)
Female Luer Thread Style Coupler



LUER