

**Material Suffix Key**

- 1 White Nylon
- 4 Green Nylon
- 6 Natural Polypropylene
- 9 Clear Polycarbonate
- 10 Acetal
- 5 Blue Nylon
- 6005 Animal-Free Polypropylene
- 9002 Radiation Stable Polycarbonate
- 2 Black Nylon
- 06 Orange Nylon
- 8003 Natural ABS
- J1A Natural Kynar PVDF
- 3 Red Nylon
- 07 Yellow Nylon
- 81 Lace White ABS

**FTLC (-1 -6 -9 -9002 -J1A)**  
Female Luer Lug Style Coupler



**FTLE (-1 -6 -9 -J1A)**  
Female Luer Lug Style Elbow



**LC34 (-1 -6 -9 -9002 -J1A)**  
Male Luer Integral Lock Ring to Female Luer Thread Style Coupler




**LE87 (-1 -6 -8003 -9 -9002)**  
Male Luer to Female Luer Thread Style Elbow (May be used with separate rotating lock ring; FSLLR)



**LC78 (-1 -6 -9 -9002 -J1A)**  
Male Luer to Female Luer Thread Style Coupler (May be used with separate rotating lock ring; FSLLR)




**LT878 (-1 -6 -9 -9002)**  
Female Luer Thread Style (2) to Male Luer Tee (May be used with separate rotating lock ring; FSLLR)



**LE7220 (-1 -2)**  
Male Luer to 200 Series Barb Elbow, 3/32" (2.4 mm) ID Tubing (May be used with separate rotating lock ring; FSLLR)




**LT787 (-1 -6 -9 -9002)**  
Male Luer (2) to Female Luer Thread Style Tee (May be used with separate rotating lock ring; FSLLR)



**LE7230 (-1 -6 -9)**  
Male Luer to 200 Series Barb Elbow 1/8" (3.2 mm) ID Tubing (May be used with separate rotating lock ring; FSLLR)



**FTLT (-1 -6 -9 -J1A)**  
Female Luer Lug Style Tee



**CA091 (-001\*)**  
Male Luer Elbow to 3/16" Classic Series Barb, (CF361-9002) with Rotating Lock Ring (FSLLR-9002)



\*Radiation Stable Polycarbonate

**FTLLE (-1 -6 -9 -J1A)**  
Female Luer Thread Style Elbow

