

**LUER TO THREAD:**

*Male*

**LUER TO THREAD:**

*Female*

**Material Suffix Key**

- |                          |                                      |                         |                    |
|--------------------------|--------------------------------------|-------------------------|--------------------|
| -1 White Nylon           | -10 Acetal                           | -2 Black Nylon          | -3 Red Nylon       |
| -4 Green Nylon           | -5 Blue Nylon                        | -06 Orange Nylon        | -07 Yellow Nylon   |
| -6 Natural Polypropylene | -6005 Animal-Free Polypropylene      | -8003 Natural ABS       | -81 Lace White ABS |
| -9 Clear Polycarbonate   | -9002 Radiation Stable Polycarbonate | -J1A Natural Kynar PVDF |                    |

**14MTLL (-1 -2 -6)**

Male Luer Integral Lock Ring with 9/16" Hex to 1/4-18 NPT Thread



**XFTL (-1 -6 -9 -J1A)**

Female Luer Lug Style with 1/4" Hex to 10-32 Special Tapered Thread



**XMLRL (-1 -6 -9 -9002)**

Male Luer with 1/4" Hex to 10-32 Special Tapered Thread (May be used with separate rotating lock ring; FSLLR)



**KFTL (-1 -6 -9 -J1A)**

Female Luer Lug Style with 1/4" Hex to 10-32 UNF Thread



**SMLRL (-1 -6 -9 -J1A)**

Male Luer with 5/16" Hex to 1/4-28 UNF Thread (May be used with separate rotating lock ring; FSLLR)



**SFTL (-1 -6 -9 -J1A)**

Female Luer Lug Style with 5/16" Hex to 1/4-28 UNF Thread



**BSMLRL (-6 -J1A)**

Male Luer with 5/16" Hex to 1/4-28 Bottom-Sealing Thread (May be used with separate rotating lock ring; FSLLR)



**18FTLL (-1 -6 -J1A)**

Female Luer Thread Style with 7/16" Hex to 1/8-27 NPT Thread



**LUER TO THREAD:**

*Female*

**18FTLL (-1 -4 -5 -6 -J1A -62\*)**

Female Luer Thread Style Elbow with 7/16" Hex to 1/8-27 NPT Thread



\*Black Polypropylene

**XFTLL (-1 -6 -9 -J1A)**

Female Luer Thread Style with 1/4" Hex to 10-32 Special Tapered Thread



**PANEL MOUNTS:**

*Female Thread Style and Lug Style Luer*

**KFTLL (-1 -6 -9 -J1A)**

Female Luer Thread Style with 1/4" Hex to 10-32 UNF Thread



**FLLB210 (-1 -2 -6 -9 -J1A)**

Female Luer Thread Style Panel Mount 1/4-28 UNF to 200 Series Barb, 1/16" (1.6 mm) ID Tubing



**SFTLL (-1 -2 -6 -9 -J1A)**

Female Luer Thread Style with 5/16" Hex to 1/4-28 UNF Thread



**FLLB220 (-1 -2 -6 -9 -J1A)**

Female Luer Thread Style Panel Mount 1/4-28 UNF to 200 Series Barb, 3/32" (2.4 mm) ID Tubing



**BSFTLL (-1 -6 -J1A)**

Female Luer Thread Style with 5/16" Hex to 1/4-28 UNF Bottom Sealing Thread



**FLLB230 (-1 -2 -6 -9 -J1A)**

Female Luer Thread Style Panel Mount 1/4-28 UNF to 200 Series Barb, 1/8" (3.2 mm) ID Tubing

