

VPSS401068N (Polycarbonate with Silicone Disk)

Check Valve, Male/Female Luer, 420 mL/min flow rate at 1 M head 30 psi back pressure, 1/4 psi inlet cracking pressure



VPSS401020N (Polycarbonate with Silicone Disk)

Check Valve, Male/Female Luer, 129 mL/min flow rate at 100mmHg, 45 psi back pressure, 1-5 psi inlet cracking pressure



VPSS401036SN (Polycarbonate with Silicone Disk)

Reflux Check Valve, Male/Female Luer, 350mL/min flow Rate at 1 M head, 45 psi back pressure, opens at 284 mmHg



VPSS401001N (ABS Housing with Silicone Disk)

Check Valve, Male Slip Luer, 420 mL/min flow rate at 1 M head, 30 psi back pressure, 1/4 psi inlet cracking pressure



VPM2930601N (White ABS)

Tube Clamp with split, 1/8" (3.2 mm) OD Tubing



VPS3901074N (Polyisoprene with ABS Acrylic Body)

IN-4000 Intermittent Injection Site to Male Luer Lock



VPM1060200N (Clear Low Density Polyethylene)

Long Spike Guard Closed at Tip



VPM0480201N (White Polyethylene)

Luer Tip Syringe Cap



VPB0850400N (Clear ABS)

Female Luer Thread Style Coupler



VPM1080303N (Blue Polyethylene)

Male Luer Integral Lock Ring Plug



VPB1000055N (Polycarbonate Body with Polyethylene Handle)

One way Stopcock, Male Slip Luer to Capped Female Luer Thread Style



VP455980 (Polycarbonate Body with Polyethylene Handle)

One Way Stopcock, Capped Male Luer Lock to Capped Female Luer Thread Style - Sterile Version



VPB1000079N (Polycarbonate Body with Polyethylene Handle)

Three Way Stopcock, (2) Capped Female Luer Thread Style and Capped Male Luer Lock



VPB15311LG2 (Rad. Stable Polycarbonate Body with White HDPE Handle)

3-Way Stopcock, (2) Female Luer Thread Style & (1) Male Luer with Rotating Lock Ring (NO caps)



VPB1000062N (Polycarbonate Body with Polyethylene Handle)

Three Way Stopcock, (2) Female Luer Thread Style and Male Slip Luer

